

DRI[®] Ethyl Alcohol Cross-Reactivity Tables Test Code # 136/136Q

POSITIVE COMPOUNDS

The following compounds tested POSITIVE using the DRI[®] DAU Ethyl Alcohol assay at the 100 mg/dL cutoff.

Positive Compounds	Trade Name	Concentration Tested (ng/mL)
Ethanol		100
n-Propanol		2,000

NEGATIVE COMPOUNDS

The following compounds tested NEGATIVE on the DRI[®] DAU Oxycodone assay at the 100 ng/mL cutoff.

Negative Compounds	Trade Name	Concentration Tested (µg/mL)
Acetaldehyde		2,000
Acetone		2,000
n-Butanol		2,000
Ethylene Glycol		2,000
Isopropanol		2,000
Methanol		2,000