



4081 Southwest 47 Street
 Suites 1-4
 Davie, FL 33314
 Tel: (954) 43 Fax: (954) 430-5715

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
% FREE PSA	%FPSA	RED	SERUM	RF
17-KETOGENIC STEROID, 24 HR URINE	4067L	24 HR UC	URINE	RF
17-KETOGENIC STEROID, URINE	4066L	UC	URINE	RF
17-KETOSTEROID, 24 HR URINE	4069L	24 HR UC	URINE	RF
17-KETOSTEROID, URINE	4068L	UC	URINE	RF
17-OH PROGESTERONE	937L	SST	SERUM	RF
17-OH/CREATININE RATIO	4065L	UC	URINE	RF
17-OH-CORTICOSTEROIDS, 24 HR URINE	986L	24 HR UC	URINE	RF
17-OH-CORTICOSTEROIDS, URINE	985L	SC	URINE	RF
17-OHPREGNENOLONE	4679L	RED SST LAV FZ	SERUM/PLASMA	FZ
5' NUCLEOTIDASE	854L	SST	SERUM	RF
5-HIAA URINE, RANDOM	4490L	UC	URINE	RF
5-HIAA, URINE 24hr	952L	24hr UC	URINE	RF
A WOMENS CENTER WEIGHT LOSS PROFILE	30781	SST,LAV,UC		
A/G RATIO	6011B	SST	SERUM	RF
A/G RATIO	A/G	SST	SERUM	RF
A/G RATIO	A/GJ	SST	SERUM	RF
AATYP DRG PNL/GY,SST	38344	GY,SST	PLASMA,BLOOD	
ABACOA PAIN PANEL	91411	SC	URINE	RF
ABE4427g MULTI DOSE STAGE, FEMALE	8911			
ABE4427g MULTI DOSE STAGE, MALE	8910			
ABE4427G-MONITORING	89101	SST,B,LAV,SC	BLD,PLASMA	
ABE4427G-SCREENING	89100	LV,SST,B,SC	BLD,PLSMA	
ABO GROUP	123L	LAV	WHOLE BLOOD	RF
ABO GROUP & RH TYPE	119B	LAV/RED	WHOLE BLOOD	RF
ABO GROUP & RH TYPE	119L	LAV/RED	WHOLE BLOOD	RF
ABSOLUTE NEUTROPHIL COUNT	246	LAV	WHOLE BLOOD	RF
ACETAMINOPHEN, SERUM	091L	SST	SERUM	RF
ACETAMINOPHEN, URINE	858M	SC	URINE	RF
ACETONE, SERUM	026L	SST	SERUM	RF
ACETONE, URINE	694L	SC	URINE	RF
ACHR BINDING Abs	580L		SERUM	RF
ACHR BLOCKING Abs	4062L	SST	SERUM	RF
ACHR MODULATING Abs	4063L	SST	SERUM	RF
ACID PHOSPHATASE, TOTAL	4289L	RED SST	SERUM	RF
ACR INITIAL PROFILE	30641	2L,1Y,1U,3SST,1	SERUM,PLAS,U	RF,FZ
ACR SURVEILLANCE	30642	SL,1Y	PLASMA	RM,FZ
ACTH, PLASMA	117L	LAV	PLASMA	RF
ACTIVATED PROTEIN C RESISTANCE	5038L	BLUE	PLASMA	FZ
ACYLCARNITINE PROFILE, PLASMA	1048L	LAV OR GREEN	PLASMA	FZ
ADENOVIRUS, IgG	705L	SST	SERUM	RF
ADENOVIRUS, IgM	706L	SST	SERUM	RF
ADH (ANTIDIURETIC HORMONE)	384L	RED	SERUM	RF
ADIPONECTIN	4290L	RED SST	SERUM	RF
ADMIT PNL-2SST/LAV	70112	SST,LAV	SERUM,PLASMA	
ADRENAL Abs, QUANT	4684L	RED SST	SERUM	RF
AFB CULTURE AND SMEAR	088L	SC-SWAB	FLUID	RF
AFB CULTURE AND SMEAR	088Q	SC	FLUID STOOL	RF
AFP Mom VALUE	456L	SST	SERUM	RF
AFP TETRA	2030L	SST	SERUM	RF
AFP VALUE (EIA)	455L	SST	SERUM	RF
AFP, MATERNAL SCREEN	468L	SST	SERUM	RF
AFP, MATERNAL SERUM	468Q	RED SST	SERUM	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
AFP, TUMOR MARKER	163	SST	SERUM	RF
AGE	491	SST	SERUM	RF
AIM PANEL-2SST-LAV	31341	SST,LAV	PLASMA,BLOOD	
ALBUMIN	124	SST	SERUM	RF
ALBUMIN	124H	SST	SERUM	RF
ALBUMIN	124J	SST	SERUM	RF
ALBUMIN	6010B	SST	SERUM	RF
ALCOHOL BLOOD	190M	RED	SERUM	RF
ALCOHOL, BLOOD	190B	SST	SST	RF
ALCOHOL, URINE (mg/dL)	136Q	UC	URINE	RF
ALDOLASE	142L	SST	SERUM	RF
ALDOSTERONE, 24HR URINE	743L	24 HR U.C.	URINE	RF
ALDOSTERONE, SERUM	143L	SST	SERUM	RM
ALDOSTERONE, URINE	741L	SC	URINE	RF
ALDOSTERONE/RENIN RATIO PNL	4708L	RED LAV SST	SERUM PLASMA	FZ
ALK PHOS	125	SST	SERUM	RF
ALK PHOS	125H	SST	SERUM	RF
ALK PHOS	125J	SST	SERUM	RF
ALK PHOS	6012B	SST	SERUM	RF
ALK PHOS ISOENZYMES	366L	SST	SERUM	RF
ALK PHOS, BONE SPECIFIC	373L	RED	SERUM	RF
ALKALINE PHOSPHATASE ISOENZYMES	366Q	RED SST	SERUM	RF
ALLAHRAKHA PANEL	3034	SST,LAV	SERUM,PLASMA	RF
ALLERGEN CLASS DESCRIPTION	2037L		SERUM	RF
ALLERGEN PANEL, BIRD	4259L	SST RED	SERUM	RF
ALLERGEN PANEL, FOOD-BASIC	4108L	SST	SERUM	RF
ALLERGEN PANEL, FOOD-FRUIT	4036L	SST	SERUM	RF
ALLERGEN PANEL, FOOD-MILK (DAIRY)	4553L	RED SST	SERUM	RF
ALLERGEN PANEL, HYMENOPTERA II	4109L	SST	SERUM	RF
ALLERGEN PANEL, MOLD	4227L	SST	SERUM	RF
ALLERGEN PANEL, SHELLFISH	3046L	SST	SERUM	RF
ALLERGEN PANEL, SPRING TREE	4354L	SST	SERUM	RF
ALLERGEN PANEL, SUMMER GRASSES	4363L	SST	SERUM	RF
ALLERGEN PANEL, VEGETABLE I	4015L	SST	SERUM	RF
ALLERGEN PANEL, VEGETABLE II	4023L	SST	SERUM	RF
ALLERGEN PANEL, ZONE 15	3026L	SST	SERUM	RF
ALLERGEN PANEL,BASIC FOOD	2049L	SST	SERUM	RF
ALLERGEN PANEL,ZONE 4	2082L	SST	SERUM	RF
ALLERGEN PNL, ZONE 15	4683L	SST	SERUM	RF
ALLERGEN PNL,ZONE 5- 2SST	2050L	2SST	SERUM	RF
ALLERGEN, INHALANTS (MIA)	4558L			
ALLERGY TEST	001L	SST	SERUM	RF
ALOPECIA PANEL	38391	2SST,LAV	PLASMA,SERUM	
ALPHA 2-MACROGLOBULINS	319L	SST	SERUM	RF
ALPHA HYDROXY ALPRAZOLAM	870M	SC	URINE	RF
ALPHA-1-ANTITRYPSIN, SERUM	160L	RED	SERUM	RF
ALPHA-GALACTOSIDASE A DEFICIENCY	4374L	ACD YELLOW	WHOLE BLOOD	RF
ALPRAZOLAM, SERUM	962M	RED GREEN	SERUM PLASMA	RF
ALPRAZOLAM, URINE	869M	SC	URINE	RF
ALT	154	SST	SERUM	RF
ALT	154H	SST	SERUM	RF
ALT	154J	SST	SERUM	RF
ALT	6014B	SST	SERUM	RF
ALT (SGPT)	154L	SST	SERUM	RF
ALTERNATIVE ADM PNL-SST,L	71201	SST LAV	blood, srm	
ALTERNATIVE DRUG SCREEN	71200	SC	URINE	RF
ALUMINUM, BLOOD	165L	ROYAL BLUE	PLASMA	RM
AMENORRHEA PROFILE	1028	SST	SERUM	RF
AMIKACIN, PEAK	186	RED	SERUM	RF
AMIKACIN, RANDOM	184	RED	SERUM	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
AMIKACIN, TROUGH	185	RED	SERUM	RF
AMINO ACID PANEL, BLOOD	1000L	GREEN	PLASMA	FZ
AMINO ACID PANEL, URINE	4136L	UC	URINE	FZ
AMINO ACID, BLOOD	1000Q	GR LAV FROZEN	PLASMA	FZ
AMINO ACID, URINE	4136Q	SC FROZEN	URINE	RF
AMINOLEVULINIC ACID, URINE	956L	SC	URINE	RF
AMIODARONE (CORDARONE), SERUM	169L	RED	SERUM	RF
AMITRIPTYLINE (ELAVIL)	848M	RED	SERUM	RF
AMITRIPTYLINE (ELAVIL)	864L	RED/LAV/GR	SERUM/PLASMA	RF
AML CLIENT PANEL-2SST*L	31720	SST,L	SERUM,PLASMA	
AML LP PANEL-2SST	31722	2SST	SERUM	
AML STD PANEL-3SST,SC	31723	3SST, SC	SERUM, URINE	
AMMONIA	172	LAV	PLASMA	FZ
AMMONIA	172B	LAV	PLASMA	FZ
AMMONIA	172J	LAV	PLASMA	FZ
AMMONIA	172L	LAV	PLASMA	FZ
AMORPHOUS	U23	SC	URINE	RF
AMPHETAMINE	802M	SC	URINE	RF
AMPHETAMINE CONF (GC/MS)	500QU		URINE	RF
AMPHETAMINE CONFIRMATION	500M	SC	URINE	RF
AMPHETAMINE CONFIRMATION, BLOOD	3079L	GRAY	WHOLE BLOOD	RF
AMPHETAMINE CONFIRMATION, URINE	3080L	SC	URINE	RF
AMPHETAMINES	500	SC	URINE	RM
AMPHETAMINES	500Q	SC	URINE	RF
AMYLASE	1250	SST	SERUM	RF
AMYLASE ISOENZYMES	4472L	RED SST	SERUM	RT
AMYLASE ISOENZYMES	4472Q	RED SST	SERUM	RF
AMYLASE, 24 HOUR URINE	4612L	24HR UC	URINE	RF
AMYLASE, BODY FLUID	4560L	SC	FLUID	RF
AMYLASE, URINE	4613L	SC	URINE	RF
ANA	045	SST	SERUM	RF
ANA CONFIRMATORY PANEL	714	2 SST	SERUM	RF
ANA DIRECT	060L	SST	SERUM	RF
ANA IFA	419L	SST	SERUM	RF
ANA SCREEN RFX TITER/PATTERN	4135Q	RED SST	SERUM	RF
ANA TITER & PATTERN	4135L	SST	SERUM	RF
ANA w/REFLEX IF POS	691	2 SST	SERUM	RF
ANA w/REFLEX IF POSITIVE	4216L	SST	SERUM	RF
ANABOLIC STEROID PNL, URINE	4464L	SC	URINE	RF
ANAEROBIC CULTURE	080	SC-SWAB	FLUID	RM
ANCA PANEL	304L	SST	SERUM	RF
ANCA PANEL	304Q	RED SST	SERUM	RF
ANDROSTANE DIOL GLUC	4685L	RED SST	SERUM	RF
ANDROSTENEDIONE	090L	SST,RED,LAV	SERUM,PLASMA	RF
ANEMIA PROFILE	2021	SST LAV	SERUM/PLASMA	RF
ANGELMAN&PRADER-WILLI SYNDROME DNA	4321L	LAV ACD	WHOLE BLOOD	RT
ANGIOTENSIN-CONVERTING ENZYME	517L	SST	SERUM	RF
ANION GAP	ANION	SST	SERUM	RF
ANISO	AMIS	LAV	WHOLE BLOOD	RF
ANTI CENTROMERE B	4218L	SST	SERUM	RF
ANTI dsDNA Ab	312	SST	SERUM	RF
ANTI GD1a, IgG	4621L	RED SST	SERUM	RF
ANTI GD1b, IgG	4623L	RED SST	SERUM	RF
ANTI GD1b, IgM	4622L	RED SST	SERUM	RF
ANTI-ASIALO GM1, IgG	4575L	RED SST	SERUM	RF
ANTI-ASIALO,GM1 IgM	4574L	RED SST	SERUM	RF
ANTIBODY ID	403B	LAV/RED	WHOLE BLOOD	RF
ANTIBODY ID #1	4370L	RED AND LAV	WB AND SERUM	RF
ANTIBODY ID #2	4595L	RED LAV	WHOLE BLOOD	RF
ANTIBODY IDENTIFICATION	4369L	RED AND LAB	WB AND SERUM	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
ANTIBODY SCREEN	407L	RED	SERUM	RF
ANTIBODY SCREEN, RFX ID	5460Q	LAV ACD	WHOLE BLOOD	RF
ANTICARDIOLIPIN Ab PROFILE	255L	RED	SERUM	RF
ANTICARDIOLIPIN IgA	254L	RED	SERUM	RF
ANTICARDIOLIPIN IgG	252L	RED	SERUM	RF
ANTICARDIOLIPIN IgM	253L	RED	SERUM	RF
ANTICENTROMERE B Ab	682L	RED	SERUM	RF
ANTICHROMATIN IgG	4329L	RED SST	SERUM	RF
ANTI-DNase(B)STREP Ab	4440L	RED SST	SERUM	RM
ANTIENDOMYSIAL Ab, IgA	678L	RED	SERUM	RF
ANTIEXTRACTABLE NUCLEAR AG	444L	SST	SERUM	RF
ANTIEXTRACTABLE NUCLEAR Ag (ENA)	371	SST	SERUM	RF
ANTI-GD1a, IgM	4620L	RED SST	SERUM	RF
ANTIGEN TYPING, RBC	184L	LAV	WHOLE BLOOD	RF
ANTIGLIADIN Abs, IgG	670L	SST	SERUM	RF
ANTIGLOMERULAR BM Ab, QUANT	317L	SST	SERUM	RF
ANTI-GM1, IgG	5028L	SST	SERUM	RF
ANTI-GM1, IgM	5029L	SST	SERUM	RF
ANTI-HISTONE Ab	315L	SST	SERUM	RF
ANTI-Hu Abs	4569L	RED SST	SERUM	RF
ANTI-INTERMYOFIBRILLAR	644L	SST	SERUM	RF
ANTIMULLERIAN HORMONE (AMH)	4319L	RED SST	SERUM	RF
ANTI-MYELOPEROXIDASE (MPO)	587L	SST	SERUM	RF
ANTIMYOCARDIAL ANTIBODIES, QN	645L	SST	SERUM	RF
ANTIPANCREATIC ISLET CELLS	4689L	RED SST	SERUM	RF
ANTIPHOSPHATIDYLERINE PNL	4331L	RED SST	SERUM	RF
ANTI-PROTEINASE 3 (PR-3)	586L	SST	SERUM	RF
ANTI-Ri Abs	4571L	RED SST	SERUM	RF
ANTI-RIBOSOMAL P ABs	309L	SST	SERUM	RF
ANTI-SARCOLEMMMA Ab	641L	SST	SERUM	RF
ANTISKIN AUTOANTIBODIES	707L	RED/SST	SERUM	RT
ANTI-SPERMATOZOA Ab	4200L	SST	SERUM	RF
ANTI-ssDNA, IgG	922L	SST RED	SERUM	RT
ANTI-STRIATION Abs	354L	SST	SERUM	RF
ANTITHROMBIN ACTIVITY	5036L	BLUE	PLASMA	FZ
ANTITHROMBIN ANTIGEN	296L	BLUE	PLASMA	FZ
ANTITHROMBIN III ACTIVITY	293L	BLUE	PLASMA	FZ
ANTITHROMBIN III PANEL	4219L	BLUE	PLASMA	FZ
ANTI-Yo Abs	4570L	RED SST	SERUM	RF
AOTRC TDM PROFILE	3032		SST,LAV	
APOLIPOPROTEIN A-1	4181L	SST	SERUM	RF
APOLIPOPROTEIN A-1	619L	SST	SERUM	RF
APOLIPOPROTEIN B	4182L	SST	SERUM	RF
APOLIPOPROTEIN B	618L	SST	SERUM	RF
APOLIPROTEIN ASSESSMENT (RATIO)	4179L	SST	SERUM	RF
APPEARANCE	4081L	SST	SERUM	RF
APTT MIXING STUDIES	2032L	BLUE	PLASMA	FZ
ARC DRUG PANEL (UA)	71220		URINE	RF
ARSENIC, BLOOD	652L	ROYAL BLUE	WHOLE BLOOD	RF
ARTHRITIS PANEL	1006	SST	SERUM	RF
ARTHUR PANEL	1611	SST,B,L,U	SERUM,PLAS,U	RF
ASHKENAZI JEWISH CARRIER PANEL	4527L	LAV AND ACD	WHOLE BLOOD	RM RF
ASIALO GM1 ANTIBODY PANEL	4572L	RED SST	SERUM	RF
ASO	046	SST	SERUM	RF
AST	153	SST	SERUM	RF
AST	153H	SST	SERUM	RF
AST	153J	SST	SERUM	RF
AST	6013B	SST	SERUM	
ATC ADMIT UDS- U	31490			
ATYP LYM	VARL	LAV	WHOLE BLOOD	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
B CELL GENE REARRANGEMENT, PCR	4052L	LAV ACD	WHOLE BLOOD	RT
B/C RATIO	B/CJ	SST	SERUM	RF
B0401008 SCREENING	8914	SST LAV UA		
B12	115	SST	SERUM	RF
BACTERIA	U15	SC	URINE	RF
BACTERIAL ANTIGENS	4632L	RED SST SC	SERUM,UR,CSF	RF
BAL HARBOUR PANEL	39180	SST,LAV, B	SERUM,BLOOD	
BARBITURATE CONF (GC/MS)	501QU	SC	URINE	RF
BARBITURATE CONFIRMATION, URINE	501L	SC	URINE	RF
BARBITURATE CONFIRMATON	501M	SC	URINE	RF
BARBITURATES	501	SC	URINE	RF
BARBITURATES	501Q	SC	URINE	RF
BARBITURATES	519L	SC	URINE	RF
BARBITURATES, BLOOD (GC/MS)	3065L	GRAY	WHOLE BLOOD	RF
BARHOUSH ANNUAL	737	2SST	SERUM	RF
BARHOUSH PRENATAL PROFILE I	16	2SST,1R,2L,U	S,P,URINE	RF
BARHOUSH PRENATAL PROFILE II	17	2SST,1R,2L,U	S,P,URINE	RF
BARTONELLA ANTIBODY PANEL	4409L	RED SST	SERUM	RF
BARTONELLA DNA by PCR	576L	SST	SERUM	RF
BASIC METABOLIC PANEL	1020	SST	SERUM	RF
BASIC METABOLIC PANEL	1020H	SST	SERUM	RF
BASIC METABOLIC PANEL	1020J	SST	SERUM	RF
BASIC METABOLIC PANEL	1020MM	SST	SERUM	RF
BASIC METABOLIC PROFILE	915B	SST	SERUM	RF
BASO ABS	8019B	LAV	WHOLE BLOOD	RF
BASO%	8014B	LAV	WHOLE BLOOD	RF
BASOPHIL	BAS#J	LAV	WB	RF
BASOPHIL	BAS%J	LAV	WB	RF
BASOPHILS	BASO%	LAV	WHOLE BLOOD	RF
BASOPHILS #	BASO	LAV	WHOLE BLOOD	RF
BASOSTIP	BAST	LAV	WHOLE BLOOD	RF
BAYVIEW PNL-2SST,LAV,SC	70471	SST,LAV,URINE	PLASMA,SERUM	
BAZILE PRENATAL PANEL	30291			
BCR-ABL1 KINASE DOMAIN	2022L	LAV OR ACD	WHOLE BLOOD	RT
BCR-ABL1 TRANSCRIPT QUAL	2024L	LAV OR ACD-A	WHOLE BLOOD	RT
BCR-ABL1 TRANSCRIPT QUANT	2023L	LAV OR ACD-A	WHOLE BLOOD	RT
BD HEALTH DRUG PNL	39521	sc	URINE	
BEHAVIOR HLTH ADM PNL-SST,L,SC	71031	SST,LAV,SC	BLD,UR,PLS	
BEHAVIORAL HEALTH UDS -SC	71030	SC	URINE	
BENZENE, BLOOD	4716L	LAV UNOPENED	WHOLE BLOOD	RF
BENZO CONFIRM,SERUM	910M	RED GREEN	SERUM PLASMA	RF
BENZODIAZAPINE CONFIRMATION	502M	SC	URINE	RF
BENZODIAZEPAM SCREEN	841M	RED	SERUM	RF
BENZODIAZEPINE	502	SC	URINE	RM
BENZODIAZEPINE	502Q	SC	URINE	RF
BENZODIAZEPINE CONF(GC/MS)	502QU	SC	URINE	RF
BENZODIAZEPINE SCREEN	839M	RED	SERUM	RF
BENZODIAZEPINES	502L	SC	URINE	RF
BENZODIAZEPINES CONFIRMATION, BLOOD	3085L	GRAY	WHOLE BLOOD	RF
BENZODIAZEPINES CONFIRMATION, URINE	3086L	SC	URINE	RF
BENZODIAZEPINES, COMP QUANT	873M	SC	URINE	RF
BENZOYLECGONINE	801M	SC	URINE	RF
BETA	4083L	SST	SERUM	RF
BETA-2 GLYCOPROTEIN I IgA	861L	SST	SERUM	FZ
BETA-2 GLYCOPROTEIN I IgG	860L	SST	SERUM	FZ
BETA-2 GLYCOPROTEIN I IgM	862L	SST	SERUM	FZ
BETA-2 MICROGLOBULIN, SERUM	186L	SST	SERUM	RF
BETA-HYDROXYBUTYRIC ACID	4427L	RED SST LAV GR	SERUM PLASMA	RF
BILE ACIDS	4439L	RED SST	SERUM	FZ
BILI	U04	SC	URINE	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
BILIRUBIN CONFIRMATION	775	SC	URINE	RF
BILIRUBIN DIRECT	128	SST	SERUM	RF
BILIRUBIN DIRECT(c)	128H	SST	SERUM	RF
BILIRUBIN INDIRECT	IBIL	SST	SERUM	RF
BILIRUBIN INDIRECT(u)	IBILH	SST	SERUM	RF
BILIRUBIN TOTAL	127	SST	SERUM	RF
BILIRUBIN TOTAL	127H	SST	SERUM	RF
BILIRUBIN, TOTAL	417L	SST	SERUM	RF
BILIRUBIN, TOTAL	6015B	SST	SERUM	RF
BIO-HEALTH PROFILE-2SST,LV,SC	32180	SST,SC,LAV	PLSM,BLD,UR	
BIOLISTIC HBV PNL-7.5 ML SRM	25101	SST.RED	SERUM	
BIOPSY, (TISSUE)	419	B. CONTAINER	FORMALIN	RM
BIOPSY, TISSUE	419Q	SC	TISSUE	RT
BIOPSY, TISSUE (2SPECIMENS)	5747Q	FORMALIN	TISSUE	RT
BIOTINIDASE DEFICIENCY	4346L	RED SST TO TT	SERUM	FZ
BK VIRUS QUANT	1041L	UC	URINE	RF
BLADDER CANCER, FISH	4647L	SC	URINE	RF
BLAST	BLAS	LAV	WHOLE BLOOD	RF
bLASTOMYCES AbS	4406L	RED SST	SERUM	RT
BLOOD	U07	SC	URINE	RF
BLOOD CULTURE	405	YELLOW	WHOLE BLOOD	RM
BLOOD CULTURE	408B	YELLOW/BOTTLE	WHOLE BLOOD	RM
BLOOM SYNDROME, DNA ANALYSIS	4545L	ACD LAV	WHOLE BLOOD	RF RM
BLUMENTHAL FEMALE PNL	39140	PROBETEC,THINP	FLUID	
BODY FLUID CULTURE	613L	SWAB	SWAB	RF
BODY FLUID CULTURE	613Q	SWAB SC	SWAB SC	RF
BODY SURFACE AREA	679	SC	URINE	RF
BORDETELLA PERTUSSIS, IgG	4337L	RED SST	SERUM	RT
BRAIN NATRIURETIC PEPTIDE	488	LAV	FROZEN PLASM	FZ
BRAIN NATRIURETIC PEPTIDE	488B	LAV	FROZEN PLASM	FR
BRAIN NATRIURETIC PEPTIDE	488L	TT	FROZEN PLASM	FR
BRIO PROFILE T3,T4,TSH,FREET4	4085	SST	SERUM	RF
BUN	129	SST	SERUM	RF
BUN	129H	SST	SERUM	RT
BUN	129J	SST	SERUM	RF
BUN	6005B	SST	SERUM	RF
BUN/CREAT RATIO	B/C	SST	SERUM	RF
BUPRENORPHINE	521N	UC		
BUPRENORPHINE	521Q	SC	URINE	RF
BUPRENORPHINE CONFIRMATION	900M	SC	URINE	RF
BUSPIRONE(BUSPAR), SERUM	970M	RED GREEN	SERUM PLASMA	RF
C&P DRUG SCREEN PANEL	88601	URINE	URINE	
C. Pneumoniae IgG	475L	SST	SERUM	RF
C. Psittaci IgG	476L	SST	SERUM	RF
C. Trachomatis IgG	477L	SST	SERUM	RF
C.A.R.E. UDS PANEL	38591	SC	URINE	
C001 PENICILLIN G	4343L	RED SST	SERUM	RF
C002 PENICILLIN V	4344L	RED SST	SERUM	RF
C1 ESTERASE INHIBITOR	814L	SST	SERUM	RF
C1 ESTERASE INHIBITOR, FUNC	4603L	RED LAV FROZEN	SERUM PLASMA	FZ
C212 ERYTHROMYCIN	3027L	SST	SERUM	RF
CA 15-3	222	SST	SERUM	RF
CA 19-9	624	SST	SERUM	RF
CA 27.29 (ICMA)	522L	SST	SERUM	RF
CA-125	316	SST	SERUM	RF
CADET PROFILE	1069	2SST	SERUM	RM
CADMIUM, BLOOD	166L	ROYAL BLUE	WHOLE BLOOD	RT
CADMIUM, URINE	4381L	SC	URINE	RF
CALCITONIN, SERUM	713L	RED TO TT	SERUM	FZ
CALCIUM	130	SST	SERUM	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
CALCIUM	130H	SST	SERUM	RF
CALCIUM	130J	SST	SERUM	RF
CALCIUM	6007B	SST	SERUM	RF
CALCIUM, 24HR URINE	625L	24 HR U.C.	URINE	RF
CALCIUM, IONIZED	204L	SST	SERUM	RF
CALCIUM, URINE	622L	SC	URINE	RF
CALL RESULT	CALL			
CANAVAN DISEASE PANEL	5009L	LAV ACD	WHOLE BLOOD	RT
CANCELLED CBC (CLOT)	1227			
CANCELLED CHEM (HEMGR)	1240			
CANCELLED COAG (CLOT)	1290			
CANCELLED COAG (HEMGR)	1243			
CANCELLED ESR (CLOT)	1230			
CANCELLED RETIC (CLOT)	1228			
CANCELLED TEST (ILLL)	1253			
CANCELLED TEST (ILLR)	1256			
CANCELLED TEST (NOICE)	1261			
CANCELLED TEST (NOID)	1249			
CANCELLED TEST (NSR)	1264			
CANCELLED TEST (NTO)	1280			
CANCELLED TEST (NTO-BLUE)	1300			
CANCELLED TEST (NTO-GR)	1310			
CANCELLED TEST (NTO-LAV)	1302			
CANCELLED TEST (NTO-RED)	1304			
CANCELLED TEST (NTO-SST)	1306			
CANCELLED TEST (NTO-UA)	1308			
CANCELLED TEST (OLD)	1267			
CANCELLED TEST (QNS)	1270			
CANCELLED TEST (UDS-NotAvail)	1315			
CANCELLED TEST (UDS-Pat Ref)	1318			
CANCELLED TEST (UDS-Unable)	1322			
CANCELLED TEST (WSR)	1275			
CANCELLED URINE (LEAK)	1246			
CANCELLED GLYCO (REFLEX)	1235			
CANDIDA ANTIBODIES, QUAL	438L	SST	SERUM	RF
CANNABINOID CONFIRMATION	510L	SC	URINE	RF
CANNABINOID CONFIRMATION, BLOOD	3090L	GRAY	WHOLE BLOOD	RF
CANNABINOID CONFIRMATION, URINE	3094L	SC	URINE	RF
CANNABINOIDS	506	SC	URINE	RF
CANNABINOIDS	506Q	SC	URINE	RF
CANNABINOIDS CONF (GC/MS)	506QU	SC	URINE	RF
CANNABINOIDS CONFIRMATION	506M	SC	URINE	RF
CANNABINOIDS, GC/MS CONFIRMATION	942M	RED GREEN	SERUM PLASMA	RF
CA-OXA	U25	SC	URINE	RF
CARBAMAZEPINE (TEGRETOL)	078	RED	SERUM	RF
CARBON DIOXIDE	126	SST	SERUM	RF
CARBON DIOXIDE	126J	SST	SERUM	RF
CARBON MONOXIDE, BLOOD	206L	ROYAL BLUE	WHOLE BLOOD	RT
CARDIOLIPIN Ab PANEL	255Q	BLUE TOP PLASMA		
CARDIOLIPIN Ab, IgA	254Q	BLUE TOP PLASMA		
CARDIOLIPIN Ab, IgG	252Q	BLUE PLASAM	PLASMA	FZ
CARDIOLIPIN Ab, IgM	253Q	BLUE TOP PLASMA		
CARISOPRODOL (SOMA) PANEL	5023L	RED LAV GREEN	SERUM PLASMA	RT
CARNITINE, TOTAL AND FREE	4301L	RED SST	SERUM	FZ
CAROTENE, BETA	209L	SST	SERUM	FZ
CAST	U22	SC	URINE	RF
CATECHOLAMINE FRACTIONATED, PLASMA	211L	GREEN	PLASMA	FROZ
CATECHOLAMINES, URINE, FREE, 24 HOUR	210L	24 HR URINE	24 HR URINE	RF
CBC	028	LAV	WB	RF
CBC	028H	LAV	WHOLE BLOOD	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
CBC	028J	LAV	WB	RF
CBC	028L	LAV	WHOLE BLOOD	RF
CBC	028MM	LAV	WHOLE BLOOD	RF
CBC (CUSTOM RANGES)	31093	LAV	WHOLE BLOOD	RF
CBC, CHEM 7, UA, HCG QUAL	916	L,SST,U	S-P-U	RM
CBCW/DIF,CMP,LIPID,TSH, T3, T4	2033	L,R,SST	SERUM/PLASMA	RM
CCP Ab, IgG	660	SST	SERUM	RF
CD4/CD8 RATIO PANEL W/O CBC	32181	ACD	WB	RT
CD4/CD8 RATIO PROFILE	435L	L/Y	WHOLE BLOOD	RM
CD4/CD8 RATIO PROFILE	435Q	LAV ACD	WHOLE BLOOD	RT
CEA	042	SST	SERUM	RF
CEFTRIAZONE (ROCEPHIN)	889L	RED	SERUM	FZ
CELIAC DISEASE COMPREHENSIVE	4186L	SST	SERUM	RF
CELIAC DISEASE HLA DQ ASSOCIATION	4610L	LAV	WHOLE BLOOD	RT
CELIAC DISEASE PANEL	646L	SST	SERUM	RF
CELIAC DISEASE PEDIATRIC PANEL	4285L	RED SST	SERUM	RF
CELL COUNT, BODY FLUID	701L	LAV,RED	FLUID	RF
CELL COUNT, FLUID	701Q	GREEN OR LAV	FLUID	RF
CELL COUNT, SYNOVIAL W/CRYSTALS	420L	GREEN LAV	FLUID	RF
CELL COUNT, SYNOVIAL, COMPLETE	420Q	SC	FLUID	RT
CELL COUNT, SYNOVIAL W/O CRYSTALS	4193L	LAV GREEN	FLUID	RF
CERULOPLASMIN	756L	SST	SERUM	RF
CHAIN OF CUSTODY SERVICE	0003			
CHALLENGES DRUG SCREEN PANEL	1037	SC	URINE	RF
CHEM26	3055	SST	SERUM	RF
CHEMISTRY PNL-SST	70781	SST	SERUM	
CHLAMYDIA & GC BY DNA PROBE	409	PROBETEC	FLUID	RF
CHLAMYDIA & GC BY DNA PROBE, URINE	411	SC	URINE	RF
CHLAMYDIA Ab PNL, IgG	4429L	RED SST	SERUM	RF
CHLAMYDIA Ab PNL, IgM	4275L	SST RED	SERUM	RF
CHLAMYDIA DNA PROBE	422L	GENPROBE	FLUID	RF
CHLAMYDIA IgG	389L	SST	SERUM	RT
CHLAMYDIA IgM	323L	RED	SERUM	RT
CHLAMYDIA PROBE	420B	GENPROBE	FLUID	RF
CHLAMYDIA PROBE, URINE	416B	SC	URINE	RF
CHLAMYDIA TRACHOMATIS	415	SC	URINE	RF
CHLAMYDIA TRACHOMATIS, DNA	611Q	THIN PREP	PAP VIAL	RT
CHLAMYDIA, NUCLEIC ACID	611L	SC,PROBE	URINE,FLUID	RF
CHLAMYDIA/GC AMPLIFICATION	601L	SC,PROBE	URINE/FLUID	RF
CHLAMYDIA/GC DNA PROBE	406L	GEN PROBE	FLUID	RF
CHLORDIAZEPOXIDE (LIBRIUM)	936M	RED GR W/FOIL	SERUM PLASMA	RF
CHLORDIAZEPOXIDE,SERUM	960L	RED	SERUM	RF
CHLORIDE	131	SST	SERUM	RF
CHLORIDE	131H	SST	SERUM	RF
CHLORIDE	131J	SST	SERUM	RF
CHLORIDE	6002B	SST	SERUM	RF
CHLORIDE, URINE	481L	SC	URINE	RF
CHLORIDE, URINE 24 HOUR	4061L	24 HR UC	URINE	RF
CHLORPROMAZINE(THORAZINE)	161M	RED	SERUM	RF
CHLORPROMAZINE, SERUM	161L	RED	SERUM	RF
CHOL / HDL	CHDL	SST	SERUM	RF
CHOLESTEROL	132H	SST	SERUM	RF
CHOLESTEROL, TOTAL	132	SST	SERUM	RF
CHOLESTEROL, TOTAL	140L	SST	SERUM	RF
CHOLESTEROL, TOTAL	4040L	SST	SERUM	RF
CHOLINESTERASE, PLASMA & RBC	2084L	LAV	WHOLE BLOOD	RF
CHOLINESTERASE, SERUM	221L	SST	SERUM	RF
CHROMIUM, PLASMA	4194L	RB	WHOLE BLOOD	RT
CHROMIUM, URINE	4248L	UC	URINE	RF
CHROMOGRANIN A	4121L	RED SST LAV GR	SERUM PLASMA	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
CHROMOSOME HIGH RES & FRAGILE X PNL	4433L	LAV GR ACD	WHOLE BLOOD	RT
CHROMOSOME MICROARRAY	4304L	GR ACD LAV	WHOLE BLOOD	RT
CHROMOSOME, BLOOD PANEL	870L	GREEN	WHOLE BLOOD	RT
CHROMOSOME, HIGH RESOLUTION	4422L	GREEN	WHOLE BLOOD	RT
CHYLOMICRONS	4082L	SST	SERUM	RF
CITALOPRAM (CELEXA)	950M	RED GREEN	SERUM PLASMA	RF
CITRIC ACID, 24 HR URINE	4498L	24 UC ACIDIFIED	URINE (ACID)	FZ
CK ISOENZYMES	74L	SST	SERUM	RF
CK ISOENZYMES	74Q	RED SST	SERUM	RF
CK TOTAL	133	SST	SERUM	RF
CK, TOTAL	133B	SST	SERUM	RF
CK-MB	092	RED SST	SERUM	RF
CLARITY	U02	SC	URINE	RF
CLASSIFICATION	4080L	SST	SERUM	RF
CLOMIPRAMINE, SERUM	4428L	RED LAV GR	SERUM PLASMA	RF
CLONAZEPAM, SERUM	837M	RED	SERUM	RF
CLONAZEPAM, URINE	836M	SC	URINE	RF
CLOSTRIDIUM DIFFICILE A&B TOXIN	141	SC	STOOL	FR
CLOSTRIDIUM DIFFICILE TOXIN A&B	413L		STOOL	FR
CLOZAPINE(CLOZARIL), SERUM	4226L	RED	SERUM	RF
CMP, URIC ACID,ALT,LIPIDS,PHOSPHORUS,LDH	756	2 SST	SERUM	RF
CMP,LIP,URICA,IRON,LDH,CBC,T4,UA,TSHGGTP	2009	2SST,LAV,UC	S/P/U	RF
CMP,LIPID,URIC A,IRON,LDH,CBC,T4,UA,GGTP	2008	2SST,LAV,UC	S/P/U	RF
CMP,LIPIDS,URIC A,IRON,CBC,LDH,T4, GGTP	2007	2SST,LAV	SERUM/PLASMA	RF
CMP,LIPIDS,URIC ACID,IRON,LDH,CBC,GGTP	2006	2SST,1LAV	SERUM/PLASMA	RF
CMP,LIPIDS,URIC ACID,IRON,LDH,GGTP	2005	2SST	SERUM	RF
CMV Ab IgG	631	SST	SERUM	RF
CMV Ab, IgM	634L	SST	SERUM	RF
CMV Ab, IgM	634Q	RED SST	SERUM	RF
CMV DNA PCR QUANT	4375L	LAV	PLASMA	FZ
CO2	126H	SST	SERUM	RF
CO2	6003B	SST	SERUM	RF
COCAINE	800M	SC	URINE	RF
COCAINE CONF (GC/MS)	503QU	SC	URINE	RF
COCAINE CONFIRMATION	503M	SC	URINE	RF
COCAINE CONFIRMATION, BLOOD	4008L	GRAY	WHOLE BLOOD	RF
COCAINE CONFIRMATION, URINE	4011L	SC	URINE	RF
COCAINE METABOLITE	503	SC	URINE	RF
COCAINE METABOLITE	503Q	SC	URINE	RF
COCAINE METABOLITE CONFIRMATION	3064L	SC	URINE	RF
COCAINE, GC/MS CONFIRMATION	947M	RED GREEN	SERUM PLASMA	RF
COCCIDIOIDES Abs	4407L	RED SST	SERUM	RT
CODEINE	816M	SC	URINE	RF
COENZYME Q10, TOTAL	5034L	GREEN TO TT	PLASMA	FZ PROTECT FROM LIGHT
COLCHICINE	4214L	SST	SERUM	RF
COLD AGGLUTININ TITER	324L	TT	SERUM	RF
COLLECTION DURATION	669	SC	URINE	RF
COLOR	U01	SC	URINE	RF
COMP DRG SCREEN	70641	SC	URINE	
COMP MET.(K+{<=3.0})	30861	SST	SERUM	RF
COMPLEMENT C1q	4513L	RED SST ON ICE	SERUM	FZ
COMPLEMENT C2	231L	RED	SERUM	FZ
COMPLEMENT C3	256	RED	SERUM	RF
COMPLEMENT C4	257	RED	SERUM	RF
COMPLEMENT C5	4039L	LAV	PLASMA	FZ
COMPLEMENT, TOTAL (CH50)	310L	SST	SERUM	RF
COMPREHENSIVE METABOLIC	1019J	SST	SERUM	RF
COMPREHENSIVE METABOLIC PANEL	1019	SST	SERUM	RF
COMPREHENSIVE METABOLIC PANEL	1019H	SST	SERUM	RF
COMPREHENSIVE METABOLIC PANEL	1019MM	SST	SERUM	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
COMPREHENSIVE METABOLIC PANEL NO GFR	10191	SST	SERUM	RF
COMPREHENSIVE METABOLIC PROFILE	1019B	SST	SERUM	RF
CONCENTRATION	132L	SC	SEMEN	RM
CONCENTRATION FOR INFECTIOUS AGENTS	994			
COOMBS DIRECT (DAT)	405B	RED LAV	SERUM/PLASMA	RF
COOMBS TITER #2	4596L	RED AND LAV	WHOLE BLOOD	RF
COOMBS, DIRECT	405L	LAV	WHOLE BLOOD	RF
COPPER, 24 HR URINE	573L	24 HR URINE	URINE	RF
COPPER, BLOOD	232L	RB LAV SST	WHOLE BLOOD	RT
COPPER, URINE	575L	SC	URINE	RF
CORTICOSTEROID BIND GLOB	4306L	RED SST	SERUM	FZ
CORTISOL	093	SST	SERUM	RF
CORTISOL BIND. GLOB WITH FREE CORTISOL	4714L	RED SST	SERUM	FZ
CORTISOL BINDING GLOBULIN	697L	SST	SERUM	FZ
CORTISOL STIMULATED	4706L	RED SST	SERUM	RF
CORTISOL, FREE, 24 HR URINE	332L	24 HR URINE	URINE	RF
CORTISOL, FREE, URINE	333L	SC	URINE	RF
COTININE QUANTITATION	847M	SC	URINE	RF
COTININE SCREEN, URINE	1047L		URINE	RF
COXSACKIE A Ab PANEL CF	715L		SERUM	RF
COXSACKIE VIRUS GROUP B Ab	939L	SST	SERUM	RF
C-PEPTIDE, SERUM	387L	SST	SERUM	RF
CREATININE	134H	SST	SERUM	RF
CREATININE	6006B	SST	SERUM	RF
CREATININE ADJUSTED VALUES PNL10	1031C			
CREATININE CLEARANCE (calc)	513	SST, 24Hr URINE	SERUM,URINE	RF
CREATININE CLEARANCE, 24 HOUR	237	SST, 24Hr URINE	SERUM,URINE	RF
CREATININE, SERUM	134	SST	SERUM	RF
CREATININE, SERUM	655	SST	SERUM	RF
CREATININE, SERUM	134J	SST	SERUM	RF
CREATININE, SERUM	139L	SST	SERUM	RF
CREATININE, URINE	53	SC	URINE	RF
CREATININE, URINE	561	SC	URINE	RF
CREATININE, URINE	235L	SC	URINE	RM
CREATININE, URINE, 24 HOUR	235	24 HR URINE	URINE	RF
CRP hs, QUANT (cardio)	174	SST	SERUM	RF
CRP, QUANTITATIVE	245	SST	SERUM	RF
CRP, QUANTITATIVE	099L	SST	SERUM	RF
CRYOGLOBULIN, QUALITATIVE	781L	RED OR SST	SERUM	RF
CRYPTOCOCCUS ANTIBODIES, QUANT	434L	SST	SERUM	RF
CRYPTOSPORIDIUM Ag	989	SC	STOOL	RF
CRYPTOSPORIDIUM/ISOSPORA SMEAR, STOOL	241L	SC	STOOL	RF
CRYSTAL EXAM, FLUID	763L	RED LAV	FLUID	RM
CRYSTALS	U16	SC	URINE	RF
CRYSTALS, FLUID	763Q	SC	FLUID	RT
CRYSTALS, SYNOVIAL FLUID	762Q	LAV GREEN	FLUID	RT
C-TELOPEPTIDE, SERUM	3063L	SST	SERUM	FZ
CULT ANY SOURCE ID	084	SWAB	FLUID	RM
CULTURE & SENSITIVITY(C/S)	082	CULTURETTE	FLUID	RM
CUSTOM PANEL 3086 (TEST)	30862			
CYCLOSPORA AND ISOSPORA EXAM	241Q	SC 10%FORMALIN	STOOL	RF
CYCLOSPORA SMEAR, STOOL	4129L	O&P PINK FORM	STOOL	RT
CYCLOSPORINE, BLOOD (FPIA)	541L	LAV/GREEN	WHOLE BLOOD	RF
CYSTATIN C, SERUM	5379L	SST	SERUM	RF
CYSTIC FIBROSIS PANEL	344L	LAV	PLASMA	RM
CYTOLOGY SMEAR MISCELLANEOUS	1046L	SLIDES	SMEAR	RT
D.A. PANEL 1	30691	SST,LAV	SERUM/WHOLE	RF
D.A. PANEL 2	30692	SST,LAV	SERUM/WHOLE	RF
D001 D PTERONYSSINUS	2051L	SST	SERUM	RF
D002 D FARINAE MITE	2052L	SST	SERUM	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
DAY 1	8912			
D-DIMER QUANTITATIVE	005L	BLUE	PLASMA	FR
DELRAY RECOVERY UDS (1UA)	71270	SC	URINE	RF
DEMETHYLDIAZEPAM	926M	RED GREEN	SERUM PLASMA	RF
DENGUE FEVER Ab PANEL	4130L	SST RED	SERUM	RF
DEOXYPYRIDINOLINE, URINE	4394L	SC	URINE	RF
DESIPRAMINE, SERUM	2003L	RED,LAV,GREEN	SERUM/PLASMA	RF
DEXTROMETHORPHAN, SERUM	958M	RED GREEN	SERUM PLASMA	RF
DEXTROMETHORPHAN, URINE	843M	SC	URINE	RF
DHEA	248L	SST	SERUM	RF
DHEA-SULFATE, SERUM	249	SST	SERUM	RF
DIA MoM	2029L	SST	SERUM	RF
DIA VALUE	2028L	SST	SERUM	RF
DIABETES PNL- GRY,LAV	39061	GREY, LAV	PLASMA,SRM	
DIABETIC PNL-2SST,L,G,SC	31032	SST,GR,L,SC	BL,PLA,UR	
DIAZEPAM (VALIUM)	928M	RED GREEN	SERUM PLASMA	RF
DIGOXIN	073	RED	SERUM	RF
DIGOXIN	073H	RED	SERUM	RT
DIGOXIN	073MM	RED	SERUM	RF
DIHYDROTESTOSTERONE	259L	TT	SERUM	FZ
DIPHENHYDRAMINE, URINE	842M	SC	URINE	RF
DIPHThERIA ANTITOXOID Ab	908L	RED SST	SERUM	RF
DISOPYRAMIDE, SERUM	175L	RED	SERUM	RF
DKOTE PNL/LV,ICELV,SST	38341	2LAV,SST	PLAMSA,SERUM	
DM PROFILE	2031	1SST,1L,1S.C.	S-P-U	RM
DOPAMINE, 24HR URINE	180L	24 HR URINE	24HR URINE	RF
DOPAMINE, URINE	178L	URINE	URINE	RF
DR SHORE ROUTINE PNL- SST,LV	31710	SST,LV	BLD,PLASMA	
DRUG ABUSE FND PNL- SST-LAV-UA	91060	SC, LAV, SST		
DRUG PANEL-ORCHID	34051	SC	URINE	RF
DRUG PNL 5, TRANSITIONS	4700	SC	URINE	RF
DRUG PNL 5,HARLEY/DAVIDSON	9044	SC	URINE	RF
DRUG SCR(10)+ EtG,urine	1036	UC	URINE	RF
DRUG SCREEN PROFESSIONAL HEALTH	3211	SC	URINE	RF
DRUG SCREEN (10)&ALC (TEST)	1072	SC	URINE	RF
DRUG SCREEN BLOOD (7)	1037L	GRAY	WHOLE BLOOD	RF
DRUG SCREEN PANEL (10), URINE	1031	SC	URINE	RF
DRUG SCREEN PANEL (5)	1030	SC	URINE	RF
DRUG SCREEN PANEL(10)QUANT,URINE	1033Q	SC	URINE	RF
DRUG SCREEN PNL(10), TEST	1071	UC	URINE	RF
DRUG SCREEN, 10 PANEL+ALC+OXY	1050	SC	URINE	RF
DRUG SCREEN, 5 PANEL + OXY	1051	SC	URINE	RF
DRUG SCREEN, BLOOD	895M	SST GREEN	SERUM PLASMA	RF
DRUG SCREEN, G&G HOLISTIC	30701	SC	URINE	RF
DRUG SCREEN, SERUM (10)	1038M	RED GREEN	SERUM	RF
DRUG SCREEN, SERUM (5)	961M	RED GREEN	SERUM PLASMA	RF
DRUG SCREEN, URINE + OXY	89501	SC	URINE	RF
DRUGS OF ABUSE SCREEN & ALCOHOL, URINE	1032	SC	URINE	RF
DRUGS OF ABUSE, SERUM (9)	1037Q	RED OR GRAY	SERUM PLASMA	RF
DSR (BY AGE)	464L	SST	SERUM	RF
DSR (SECOND TRIMESTER)	463L	SST	SERUM	RF
DTX UNIT ADMIT-SST*LAV	3086	SST LAV		
DUMENY PROFILE	3010	SST,S.C.	SERUM,URINE	RM
E001 CAT HAIR DANDER STD	2053L	SST	SERUM	RF
E002 DOG EPITHELIAL	2054L	SST	SERUM	RF
E003 HORSE EPITHELIAL/DANDER	4238L	SST	SERUM	RF
E007 PIGEON DROPPINGS	4257L	RED SST	SERUM	RF
E026 PARROT SERUM	4256L	RED SST	SERUM	RF
E070 GOOSE FEATHERS	4239L	SST	SERUM	RF
E077 BUDERIGAR DROPPINGS	4252L	RED SST	SERUM	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
E078 BUDERIGAR FEATHER	4253L	RED SST	SERUM	RF
E079 BUDERIGAR SERUM	4254L	RED SST	SERUM	RF
E213 PARROT FEATHERS	4255L	RED SST	SERUM	RF
EAR CULTURE	615	SWAB	SWAB	RM
EARLY TERMINATION/EARLY DISCONTINUATION	8913			
EBV DNA QUANT, PCR	4693L	LAV	PLASMA	FZ
EBV EA Ab, IgG	318L	SST	SERUM	RF
EBV EA ANTIBODY PANEL	4563L	RED SST	SERUM	RF
EBV NA Ab, IgG	425	SST	SERUM	RF
EBV VCA Ab, IgG	382	SST	SERUM	RF
EBV VCA Ab, IgM	380	SST	SERUM	RF
eCREATININE CLEARANCE	5131	RED SST	SERUM	RF
ECTASY	511	SC	URINE	RF
ECTASY	511Q	SC	URINE	RF
eGFR	495H	SST	SERUM	RF
eGFR	GFRJ	SST	SERUM	RF
eGFR(AFRICAN AMERICAN)	GFRAA	SST	SERUM	RF
eGFR(FEMALE)	GFRF	SST	SERUM	RF
eGFR(MALE)	GFR	SST	SERUM	RF
EHEC SHIGA TOXINS	592			
EHRlichia CHAFFEENIS-HME	4486L	RED SST	SERUM	RF
ELECTROLYTE PANEL	1010	SST	SERUM	RF
ELECTRONICALLY SIGNED OUT	0045			
ELECTRONICALLY SIGNED OUT	0047			
ENTAMOEBA HISTOLYTICA AG	4710L	OVP	STOOL	RF
ENVIROMENTAL CULTURE ANY ID.	3020	SWAB	SWAB	RF
EOS ABS	8018B	LAV	WHOLE BLOOD	RF
EOS%	8013B	LAV	WHOLE BLOOD	RF
EOSINOPHIL #	EOSI	LAV	WHOLE BLOOD	RF
EOSINOPHIL, URINE	710L	SC	URINE	RF
EOSINOPHILS	EOSI%	LAV	WHOLE BLOOD	RF
EOSONOPHIL COUNT, URINE	516			
EPINEPHRINE, 24HR URINE	182L	24 HR URINE	24 HR URINE	RF
EPINEPHRINE, URINE	181L	SC	URINE	RF
EPITHELIAL CELLS	U17	SC	URINE	RM
ERYTHROPOIETIN	264	SST	SERUM	RF
ESCITALOPRAM (LEXAPRO)	931M	RED GREEN	SERUM PLASMA	RF
estimated GFR	494			
estimated GFR	495	SST	SERUM	RF
estimated GFR	8681		SERUM	RF
ESTRADIOL	094	SST	SERUM	RF
ESTRADIOL, SENSITIVE	0994L	RED SST	SERUM	RF
ESTRIOL, SERUM	224L	SST	SERUM	RF
ESTROGENS, TOTAL	266L	SST	SERUM	RF
ESTROGENS, TOTAL	266Q	RED SST	SERUM	RF
ESTRONE, SERUM	268L	SST	SERUM	RF
ETHANOL, BLOOD	271L	GRAY OR SST	WB OR SERUM	RF
ETHOSUXIMIDE, SERUM	4051L	RED SST	SERUM	RF
ETHYL GLUCURONIDE (EtG)	533Q	UC	URINE	RF
EX71 FEATHER MIX	4478L	SST	SERUM	RF
EYE CULTURE	614	SWAB	SWAB	RM
F.C.1L,2SST,UA	90941	2SST,1LAV,1UA		
F001 EGG WHITE	2097L	SST	SERUM	RF
F002 MILK (COW)	2039L	SST	SERUM	RF
F003 CODFISH	3048L	SST	SERUM	RF
F004 WHEAT	2040L	SST	SERUM	RF
F008 CORN	2041L	SST	SERUM	RF
F012 PEA, GREEN	4024L	SST	SERUM	RF
F013 PEANUT	2042L	SST	SERUM	RF
F014 SOYBEAN	2043L	SST	SERUM	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
F018 BRAZIL NUT	4198L	SST	SERUM	RF
F020 ALMOND	4196L	SST	SERUM	RF
F023 CRAB	3050L	SST	SERUM	RF
F024 SHRIMP	3051L	SST	SERUM	RF
F025 TOMATO	4025L	SST	SERUM	RF
F026 PORK	2044L	SST	SERUM	RF
F027 BEEF	2045L	SST	SERUM	RF
F029 WATERMELON	4103L	SST	SERUM	RF
F031 CARROT	4026L	SST	SERUM	RF
F033 ORANGE	4618L	RED SST	SERUM	RF
F033 SCALLOP	3052L	SST	SERUM	RF
F035 POTATO, WHITE	4027L	SST	SERUM	RF
F037 MUSSEL	4690L	RED SST	SERUM	RF
F044 STRAWBERRY	4104L	SST	SERUM	RF
F046 YEAST, BAKER'S	4619L	RED SST	SERUM	RF
F048 ONION	4029L	SST	SERUM	RF
F049 APPLE	4033L	SST	SERUM	RF
F052 CHOCOLATE / COCOA	2048L	SST	SERUM	RF
F075 EGG (YOLK)	2098L	SST	SERUM	RF
F076 LACTOALBUMIN, ALPHA	4554L	RED SST	SERUM	RF
F077 LACTOGLOBULIN, BETA	4555L	RED SST	SERUM	RF
F078 CASEIN	2099L	SST	SERUM	RF
F079 GLUTEN	4397L	RED SST	SERUM	RF
F080 LOBSTER	3054L	SST	SERUM	RF
F081 CHEESE, CHEDDAR	4556L	SST RED	SERUM	RF
F082 CHEESE, BLUE (MOLD)	4557L	RED SST	SERUM	RF
F083 CHICKEN	3000L	SST	SERUM	RF
F084 KIWI FRUIT	3070L	SST	SERUM	RF
F085 CELERY	4018L	SST	SERUM	RF
F087 MELON, HONEYDEW	4105L	SST	SERUM	RF
F088 LAMB	4367L	SST	SERUM	RF
F092 BANANA	3068L	SST	SERUM	RF
F094 PEAR	4035L	SST	SERUM	RF
F095 PEACH	4034L	SST	SERUM	RF
F096 AVOCADO	3067L	SST	SERUM	RF
F201 PECAN NUT	4197L	SST	SERUM	RF
F202 CASHEW NUT	4195L	SST	SERUM	RF
F207 CLAM	3049L	SST	SERUM	RF
F212 MUSHROOM	4037L	SST	SERUM	RF
F214 SPINACH	4016L	SST	SERUM	RF
F215 LETTUCE	4019L	SST	SERUM	RF
F216 CABBAGE	4021L	SST	SERUM	RF
F225 SUMMER SQUASH	4031L	SST	SERUM	RF
F244 CUCUMBER	4022L	SST	SERUM	RF
F245 EGG, WHOLE	2047L	SST	SERUM	RF
F256 WALNUT, FOOD	4199L	SST	SERUM	RF
F259 GRAPE	4032L	SST	SERUM	RF
F260 BROCCOLI	4020L	SST	SERUM	RF
F287 BEAN, KIDNEY (RED)	4030L	SST	SERUM	RF
F290 OYSTER	3053L	SST	SERUM	RF
F291 CAULIFLOWER	4017L	SST	SERUM	RF
F299 SWEET CHESTNUT	3069L	SST	SERUM	RF
F315 BEAN, GREEN	4028L	SST	SERUM	RF
FACTOR II ACTIVITY	183L	TT	BLUE PLASMA	FZ
FACTOR II, DNA ANALYSIS	179L	LAV/ACD/SWAB	WHOLE BLOOD	RT
FACTOR IX ANTIGEN	4542L	BLUE PLASMA	PLASMA	FZ
FACTOR V ACTIVITY	658L	BLUE	PLASMA	FZ
FACTOR V LEIDEN MUTATION	365L	LAV/ACD	WHOLE BLOOD	RM
FACTOR VII ACTIVITY	115L	TT	PLASMA	FR
FACTOR VII ANTIGEN	4543L	BLUE PLASMA	PLASMA	FZ

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
FACTOR VIII ACTIVITY	189L	BLUE	PLASMA	FZ
FACTOR X ANTIGEN	4544L	BLUE PLASMA	PLASMA	FZ
FACTOR XI ACTIVITY	321L	TT	PLASMA (BLUE	FZ
FACTOR XII ACTIVITY	4602L	CITRATED PLASMA	PLASMA	FZ
FAMILIAL DYSAUTONOMIA	4415L	LAV ACD	WHOLE BLOOD	RT
FANCONI ANEMIA, TYPE C, DNA ANALYSIS	4547L	LAV ACD	WHOLE BLOOD	RF RT
FASTING PANEL (2SST)	91512			
FATTY ACIDS, FREE (NONESTERASE)	275L	SST	SERUM	RF
FAX	FAX			
FDP, PLASMA	4246L	BLUE	PLASMA	FZ
FEBRILE AGGLUTININ PANEL	4719L	RED SST	SERUM	RF
FECAL FAT, QUALITATIVE	671L	SC	STOOL	
FECAL REDUCING SUBSTANCES	672L	SC	STOOL	
Female admit panel- 2 sst + lav	30892	2 sst 1 lav		
FEMALE ADMIT PANEL- 2SST + 1LAV	1060A	SST	SERUM	RF
FEMALE ADM-SST,L,SC	70352	SST,LAV,SC	SERUM,PLASM	RM
FEMALE HORMONE-2SST	95402	SST	SERUM	
FEMALE PANEL - SST+LAV+SC	30912	SST,LAV,SC		
FEMALE PANEL- 4SST,LV,B	71042	LAV,4SST,BL	ALL	
FENTANYL, URINE	714L	SC	URINE	RT
FENTANYL, URINE	714M	SC	URINE	RF
FERRITIN	095	SST	SERUM	RF
FERRITIN(CUSTOM)	90810	SST	SERUM	RF
FETAL FIBRONECTIN	277L	SWAB	SWAB	FZ
FIBRINOGEN ACTIVITY	4209L	BLUE	PLASMA	FZ
FIBRINOGEN ANTIGEN	295L	BLUE	PLASMA	FZ
FIBRINOGEN EVALUATION PANEL	4208L	BLUE	PLASMA	FZ
FL RESEARCH CHEM PNL	39071		SERUM	
FL RESEARCH URINE PNL	39072		URINE	
FLECAINIDE (TAMBOCOR)	703L	RED	SERUM	RF
FLOW CYTOMETRY	4680L	LAV ACD	WHOLE BLOOD	RT
FLUID CYTOLOGY	409L	RED	FLUID	RF
FLUORESCENCE IN SITU HYBRIDIZATION,FISH	4057L	GREEN	WHOLE BLOOD	RT
FLUORESCENT SMEAR REPORT	942L	SWAB	SWAB	RM
FLUORIDE	278L	RED OR RB	SERUM/PLASMA	RF
FLUOXEINE (PROZAC),SERUM	972M	RED GREEN	SERUM PLASMA	RF
FLUOXETINE (PROZAC), SERUM	3073L	RED,LAV,GREEN	SERUM/PLASMA	RF
FLUPHENAZINE (PROLIXIN)	797L	RED	SERUM	RF
FLUVOXAMINE(LUVOX),SERUM	930M	RED GR	SERUM PLASMA	RF
FMALE >40-2SST,LAV,SC,HCARDS	31034			
FOLATE	096	SST	SERUM	RF
FOLATE, RBC	4383L	2 LAV	WHOLE BLOOD	1 FROZEN & 1 NON-FROZEN
FOLLOW-UP PROFILE	2202	SST,L,GY	SERUM/PLASMA	RM
FRAGILE X DNA	1090L	LAV	WB	RF
FREE KAPPA & LAMBDA LIGHT CHAINS, URINE	749L	24 HR UC	24 HR URINE	RT
FREE KAPPA & LAMBDA LIGHT CHANINS, SERUM	4389L	RED SST TO TT	SERUM TO TT	RF
FRUCTOSAMINE	331	RED SST	SERUM	RF
FRUCTOSAMINE	331L	SST,LAV	SERUM/PLASMA	RF
FSH	097	SST	SERUM	RF
FSH	097L	SST	SERUM	RF
FTA-ABS (T.PALLIDUM Ab)	401L	SST	SERUM	RF
FTA-Abs, IgG	483L	SST	SERUM	RF
FTA-Abs, IgM	484L	SST	SERUM	RF
FTI	FTI	SST	SERUM	RF
FTI	FTI2	SST	SERUM	RF
FUNGAL ANTIBODIES QUANT	4401L	RED SST	SERUM	RT
FUNGAL CULTURE	408L	SWAB	SWAB	RM
FUNGUS CULTURE, BLOOD	410Q	SC	BLOOD	RT
FUNGUS CULTURE, OTHER	408Q	SWAB, SC		

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
FUNGUS CULTURE,HAIR,SKIN,NAIL	409Q	HAIR SKIN NAIL		RT
FX02 FISH/SHELL MIX	2046L	SST	SERUM	RF
G and G Venipuncture	30700			
G&G FEMALE ADM-4SST*LAV*SC	30704	4SST,LAV,SC	SERUM,PLS,U	
G&G MALE ADM-4SST*LAV*SC	30702	2SST	SERUM	RF
G002 BERMUDA GRASS	2055L	SST	SERUM	RF
G004 FESCUE, MEADOW	4364L	SST	SERUM	RF
G005 RYE, PERENNIAL	2056L	SST	SERUM	RF
G006 TIMOTHY	4365L	SST	SERUM	RF
G008 BLUEGRASS, KENTUCKY	3016L	SST	SERUM	RF
G010 JOHNSON GRASS	2057L	SST	SERUM	RF
G017 BAHIA GRASS	2058L	SST	SERUM	RF
G-6-PD, QUANT, BLOOD&RBC	1099L	LAV	WHOLE BLOOD	RF
GABAPENTIN, SERUM	662L	RED	SERUM	RF
GAD-65 AUTOANTIBODY	4270L	RED SST	SERUM	RF
GASTRIN, SERUM	326L	RED SST	SERUM	FZ
GC DNA PROBE	421B	PROBETEC	FLUID	RF
GC DNA PROBE	423L	GENPROBE	FLUID	RF
GC PROBE, URINE	415B	SC	URINE	RF
GC/CHLAMYDIA DNA PROBE, PAP VIAL	601Q	THIN PREP	THIN PREP	RT
GENERAL HEALTH PROFILE	1014	1SST,1LAV	SERUM/PLASMA	RM
GENERAL OXIDANT	537	UC	URINE	RF
GENESIS DRUG PANEL	90881	SC	URINE	
GENETIC EXAM	988			
GENITAL CULTURE	059	SWAB	FLUID	RM
GENTAMICIN, PEAK	015	RED	SERUM	RF
GENTAMICIN, PEAK	015L	RED	SERUM	RF
GENTAMICIN, RANDOM	016	RED	SERUM	RF
GENTAMICIN, TROUGH	014	RED	SERUM	RF
GENTAMICIN, TROUGH	014L	RED	SERUM	RF
GESTAT. AGE BASED ON	449L	SST	SERUM	RF
GESTATIONAL AGE	448L	SST	SERUM	RF
GGT	137	SST	SERUM	RF
GGT	138L	SST	SERUM	RF
GIARDIA Ag	990	SC	STOOL	RF
GIARDIA AG, STOOL	361Q	STOOL FORMALIN	STOOL	RF
GIARDIA LAMBLIA Ag, EIA	361L	SC	STOOL	RM
GLIADIN Ab PANEL	4221L	SST	SERUM	RF
GLOBULIN	GLOB	SST	SERUM	RF
GLOBULIN	GLOBJ	SST	SERUM	RF
GLUCAGON, PLASMA	4235L	LAV	PLASMA	FZ
GLUCOSE	138	SST	SERUM	RF
GLUCOSE	138H	SST	SERUM	RF
GLUCOSE	138J	SST	SERUM	RF
GLUCOSE	6004B	SST	SERUM	RF
GLUCOSE	U03	SST	SERUM	RF
GLUCOSE 2H PP	139	GRAY	PLASMA	RF
GLUCOSE TOLERANCE, 2 HOUR	606	4 GRAY	PLASMA	RF
GLUCOSE TOLERANCE, 3 HOUR	607	5 GRAY	PLASMA	RF
GLUCOSE TOLERANCE, 4 HOUR	608	6 GRAY	PLASMA	RF
GLUCOSE TOLERANCE, 5 HOUR	609	7 GRAY	SERUM	RF
GLUCOSE, 1 HOUR	601	3 GRAY	PLASMA	RF
GLUCOSE, 1/2 HOUR	600	2 GRAY	PLASMA	RF
GLUCOSE, 2 HOUR	602	3 GRAY	PLASMA	RF
GLUCOSE, 3 HOUR	603	4 GRAY	PLASMA	RF
GLUCOSE, 4 HOUR	604	5 GRAY	PLASMA	RF
GLUCOSE, 5 HOUR	605	6 GRAY	PLASMA	RF
GLUCOSE, BODY FLUID	4507L	SC	BODY FLUID	RT
GLUCOSE, FASTING	599	7 GRAY	PLASMA	RF
GLUTATHIONE, BLOOD	708L	RED	SERUM	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
GLYCOHEMOGLOBIN	101L	LAV	WHOLE BLOOD	RF
GLYCOHEMOGLOBIN (HbA1c)	101	LAV	WHOLE BLOOD	RF
GLYCOMARK	4222L	SST LAV	SERUM PLASMA	RT
GLYCOSAMINOGLYCANS, URINE	1081L	SC	URINE	FZ
GM-1 ANTIBODY PANEL	5030L	SST	SERUM	RF
GONOCOCCUS, NUCLEIC ACID	612L	URINE/SWAB	URINE/SWAB	RF
gp120 Ab	5225L	SST	SERUM	RF
gp160 Ab	5224L	SST	SERUM	RF
gp40 Ab	5229L	SST	SERUM	RF
gp41 Ab	5228L	SST	SERUM	RF
GRAM STAIN	085	SLIDE	FLUID	RM
GRANULAR CAST	U19	UC	URINE	RF
GRASSI PROFILE #1 - SST(2)	90921	2 SST	PLASMA	
GRASSI PROFILE #2 - SST(2),LAV,SC	90922	SST,LAV,SC		
GROWTH HORMONE, 24HR URINE	4670L	24HR UC	URINE	RF
GROWTH HORMONE, SERUM	776L	SST	SERUM	RF
H PYLORI BREATH TEST	4691L	BREATH BAGS	BREATH	RT
H. INFLUENZAE B IgG	929L	RED SST	SERUM	RF
H. PYLORI Ab, IgG	199	SST	SERUM	RF
H. PYLORI AG, STOOL	446Q	SC	STOOL	RF
H. PYLORI IgG, Abs	199L	SST	SERUM	RF
H. PYLORI STOOL ANTIGEN	446L	SC	STOOL	FZ
H. PYLORI, IgA	684L	SST	SERUM	RF
H. PYLORI, IgG	685L	SST	SERUM	RF
H. PYLORI, IgM	686L	SST	SERUM	RF
H002 HOUSE DUST	4240L	SST	SERUM	RF
HALE PAIN PANEL	39041	SC	URINE	
HALLWAY DRUG PNL- SC	25110	SC	URINE	
HALOPERIDOL (HALDOL), SERUM	311L	RED	SERUM	
HANLEY FEMALE UDS ADMIT	98003	SC	URINE	RF
HANLEY MALE UDS ADMIT	98002	SC	URINE	RF
HANLEY SERUM/PLASMA ADMIT PANEL	98001	SST,LAV	SERUM,PLASMA	RF
HANNA COMP CYTO PNL(PAP+HPV)	39382	THIN PREP	THIN PREP	RF
HANNA CYTO PNL (PAP+HPV+GC/CHLAM)	39383	THIN PREP	THIN PREP	
HANNA OB PNL-SST,R,2LAV	39381	SST,R,LAV	SST,PLS,BL	
HAPTOGLOBIN QUANT	242L	SST	SERUM	RF
HBV DNA, QUALITATIVE PCR	765L	SST,ACD,LAV	SERUM/PLASMA	FZ
HBV DNA, QUANT RFX DRUG RESISTANCE	4269L	SST RED LAV	SERUM PLASMA	FZ
HBV DNA, QUANTITATIVE	873L	RED SST LAV	SERUM PLASMA	FRZ- THIS USE TO BE CALL HBV QUANTASURE
HBV GENOTYPE	4260L	SST LAV ACD	SERUM PLASMA	FZ
HBV SUPERQUANT (NGI)	4479L	LAV ACD PPT	PLASMA	FZ
hCG Mom	458L	SST	SERUM	RF
HCG QUALITATIVE SERUM	048L	SST	SERUM	RF
hCG VALUE	457L	SST	SERUM	RF
HCG, QUALITATIVE (POS=CRIT)	0481			
HCG, QUALITATIVE SERUM	048MM	SST RED	SERUM	RF
HCG, QUALITATIVE(SERUM)	048	SST	SERUM	RF
HCG, QUALITATIVE(SERUM)	048H	SST	SERUM	RF
HCG, QUALITATIVE(URINE)	623	SC	URINE	RF
HCG, QUALITATIVE(URINE)	623H	SC	URINE	RF
HCG, QUALITATIVE, URINE	623MM	SC	URINE	RF
HCG, QUANTITATIVE	049	SST	SERUM	RF
HCG, QUANTITATIVE, SERUM	049L	SST	SERUM	RF
HCV AB, RIBA	4617L	RED SST	SERUM	RF
HCV FIBROSURE	432L	SST	SERUM	RF
HCV GENOTYPING	2612L	SST,LAV	SERUM,PLASMA	FZ
HCV QUANTASURE	877L	SST,LAV	SERUM,PLASMA	FZ
HCV RIBA	4617Q	SST RED	SERUM	RF
HCV RNA, QUANT RT PCR	877Q	RED SST	SERUM	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
HCV SUPERQUANT (NGI)	4470L	PPT	PLASMA	FZ
HDL	144	SST	SERUM	RF
HDL	144H	SST	SERUM	RF
HDL CHOLESTEROL	4042L	SST	SERUM	RF
HDL CHOLESTEROL	4076L	SST	SERUM	RF
HE PROTEIN 4	4658L	RED SST	SERUM	RF
HEAVY METAL PROFILE I, BLOOD	664L	LAV/ROYAL BLUE	WHOLE BLOOD	RM
HEAVY METALS PNL I, URINE	663L	SC	URINE	RF
HEAVY METALS, PNL II, URINE	4382L	SC	URINE	RF
HELLER ANEMIA PROFILE	38012	SST LAV		
HELLER THYROID PANEL	38013	SST	SERUM	RF
HEMATOCRIT	287	LAV	WHOLE BLOOD	RF
HEMATOCRIT	31092	LAV	WHOLE BLOOD	RF
HEMATOCRIT	287J	LAV	WB	RF
HEMATOCRIT	287MM	LAV	WHOLE BLOOD	RF
HEMATOCRIT	8003B	LAV	WHOLE BLOOD	RF
HEMOGLOBIN	288	LAV	WHOLE BLOOD	RF
HEMOGLOBIN	31091	LAV	WHOLE BLOOD	RT
HEMOGLOBIN	288J	LAV	WB	RF
HEMOGLOBIN	288MM	LAV	WHOLE BLOOD	RF
HEMOGLOBIN	8002B	LAV	WHOLE BLOOD	RF
HEMOGLOBIN & HEMATOCRIT	135	LAV	WHOLE BLOOD	RF
HEMOGLOBIN A	80L	LAV	WHOLE BLOOD	RF
HEMOGLOBIN A2	82L	LAV	WHOLE BLOOD	RF
HEMOGLOBIN C	84L	LAV	WHOLE BLOOD	RF
HEMOGLOBIN ELECTROPHORESIS	100L	LAV	WHOLE BLOOD	RF
HEMOGLOBIN ELECTROPHORESIS	100Q	LAV	WHOLE BLOOD	RF
HEMOGLOBIN F	81L	LAV	WHOLE BLOOD	RF
HEMOGLOBIN S	83L	LAV	WHOLE BLOOD	RF
HEMOGLOBIN SOLUBILITY	85L	LAV	WHOLE BLOOD	RF
HEMOGLOBIN VARIANTS	4205L	LAV	WHOLE BLOOD	RF
HEMOGRAM	701	LAV	WHOLE BLOOD	RF
HENDERSON ADMIT PANEL	3005	SST,R,LAV	SERUM,PLASMA	RM
HEP A Ab, IgM	062	SST	SERUM	RF
HEP A Ab, TOTAL	072L	SST	SERUM	RF
HEP B c Ab, IgM	064	SST	SERUM	RF
HEP B c Ab, TOTAL	393L	SST	SERUM	RF
HEP B c IgM	392L	SST	SERUM	RF
HEP B s Ab	394	SST	SERUM	RF
HEP B s Ag	069	SST	SERUM	RF
HEP B s Ag	069L	SST	SERUM	RF
HEP B s Ag CONFIRMATION	5222L	SST	SERUM	RF
HEP B VIRUS PROFILE	509L	SST	SERUM	RF
HEP Be Ab	874L	SST	SERUM	RF
HEP Be Ag	875L	SST	SERUM	RF
HEP C Ab w/REFLEX TO RIBA	263	RED SST	SERUM	RF
HEP C Ab, IgG	065	SST	SERUM	RF
HEP C RNA-PCR, QUAL	551L	SST,LAV	SERUM,PLASMA	FZ
HEP C VIRUS Ab	5042L		SERUM	RF
HEPARIN Anti-Xa	888L	BLUE	PLASMA	FZ
HEPARIN INDUCED PLAT Ab	4279L	RED	SERUM	FZ
HEPATIC FUNCTION PANEL	1013	SST	SERUM	RF
HEPATIC FUNCTION PANEL	1013H	SST	SERUM	RF
HEPATIC FUNCTION PANEL	1013MM	SST	SERUM	RF
HEPATITIS DELTA Ab	4534L	RED SST	SERUM	FZ
HEPATITIS PANEL	1015	SST	SERUM	RF
HEREDITARY HEMOCHROMATOSIS, DNA	752L	LAV, ACD, SWAB	WB OR SWAB	RF/RT
HEROIN METABOLITE	520	SC	URINE	RF
HEXAGONAL PHASE PHOSPHOLIPID	994L	BLUE	PLASMA	FZ
HHV-6, IgG, QUANT	4417L	RED SST	SERUM	RT

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
HHV-6, IgM, QUANT	4431L	SST	SERUM	RF
HILLS ANEMIA PNL IV	30633	2SST -LAV		
HILLS ARTHRITIS PNL	30631	R,L,SST	SERUM,PLASMA	RM
HILLS CARDIAC PROFILE	30632	SST	SERUM	RF
HILLS GENERAL HLTH	30634	2SST- LAV		
HIRSUTISM PROFILE	260L	2SST	SERUM	RF
HISTAMINE, 24 HOUR URINE	4509L	24 HR UC	URINE	RF
HISTAMINE, PLASMA	4247L	LAV	PLASMA	FZ
HISTOPLASMA Abs	4408L	RED SST	SERUM	RT
HISTOPLASMA ANTIGEN	4540L	RED SST	SERUM	RF
HIV 1/0/2 Abs, QUAL	4258L	SST RED	SERUM	RF
HIV 1/2 SCREEN, EIA	071Q	SST RED	SERUM	RF
HIV Ab, HIV-1 WESTERN BLOT	5223Q	RED SST	SERUM	RF
HIV Abs/WESTERN BLOT	5223	SST	SERUM	RF
HIV DNA PCR (PROVIRAL)	4342L	ACD LAV	WHOLE BLOOD	RT
HIV GENOTYPE (GenoSURE)	2611L	2LAV,PPT,ACD	PLASMA	FROZEN
HIV PHENOTYPE	2613L	2 LAV	PLASMA	FZ
HIV-1 / HIV-2 Ab	071	SST	SERUM	RF
HIV-1 DNA PCR, QUAL	4633L	RED SST	SERUM	RF
HIV-1 RNA b-DNA (BRANCHED)	4338L	PPT ACD LAV	PLASMA	FZ
HIV-1 RNA by PCR (VIRAL LOAD)	506L	2L	PLASMA	FR
HIV-1 RNA by PCR (VIRAL LOAD)	5750Q	LAV PLASMA FZ	PLASMA	FZ
HIV-1 RT-PCR, QUANT	709L	LAV,PPT, ACD	PLASMA	FZ
HLA A,B DISEASE ASSOCIATION	4604L	LAV WHOLE BLOOD	WHOLE BLOOD	RT
HLA-B27	523L	LAV	PLASMA	RF
HLA-B5701	4286L	LAV BUCCAL SWAB	WB SWAB	RT
HMG Profile	38011	SST	SST	
HOMOCYSTEINE	429	SST,RED,GR,LAV	SERUM PLASMA	RF
HPV	916L	THP	FLUID	RM
HPV DNA HIGH RISK	914Q	THIN PREP	PAP VIAL	RT
HPV DNA HIGH/LOW	915Q	THIN PREP	PAP VIAL	RT
HPV, HIGH & LOW RISK	915X	THP	FLUID	RM
HPV, HIGH RISK	914X	THP	FLUID	RF
HPV, LOW RISK	916X	THP	FLUID	RT
HPV, LOW RISK	C5	THP	FLUID	RM
HSV CULTURE & TYPING	4585L	VTM	FLUID/SWAB	RFF
HSV CULTURE WITHOUT TYPING	412L	VTM	FLUID	RM
HSV CULTURE, RAPID	412Q	SC VTM	FLUID	RF
HSV I&II, IgM, TYPE SPECIFIC	4650L	SST RED	SERUM	RF
HSV I/II IgG combination	105L	SST	SERUM	RF
HSV I/II IgM Abs, INDIRECT	4702L	SST RED	SERUM	RF
HSV I/II IgM combination	104L	SST	SERUM	RF
HSV I/II, IgM	104Q			
HSV PROFILE	1002	SST	SERUM	RF
HSV TYPE I SPECIFIC, IgG	954	RED SST	SERUM	RF
HSV TYPE II SPECIFIC IgG	969	RED SST	SERUM	RF
HTLV-I/II Ab, QUALITATIVE	745L	RED	SERUM	RT
HUMAN ANTI-MOUSE Abs	4646L	RED SST	SERUM	RF
HUNTINGTON DISEASE PROFILE	890L			
HYALINE CAST	U18	SC	URINE	RF
HYDROCODONE, URINE	959M	UC	URINE	RF
HYDROMORPHONE, URINE	868M	SC	URINE	RF
HYDROXYPROLINE, TOTAL QN/24HOUR	251L	24 HR URINE	URINE	RF
HYDROXYZINE(VISTARIL),SERUM	920M	RED GREEN	SERUM PLASMA	RF
HYPCHROM	HYPC	LAV	WHOLE BLOOD	RF
HYPER COAGULATION PNL	2500	SST BLUE		
HYPERSENSITIVITY PNEUMONITIS	885L	SST	SERUM	RF
I001 HONEYBEE	4110L	SST	SERUM	RF
I002 HORNET, YELLOW FACE	4111L	SST	SERUM	RF
I003 YELLOW JACKET	4112L	SST	SERUM	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
I004 WASP, PAPER	4113L	SST	SERUM	RF
I005 HORNET, YELLOW	4114L	SST	SERUM	RF
I054 FLEA (CAT / DOG)	2083L	SST	SERUM	RF
I080 BUMBLEBEE	4115L	SST	SERUM	RF
I100 COCKROACH, AMERICAN	2059L	SST	SERUM	RF
IA-2 AUTOANTIBODIES	4426L	RED SST GR LAV	SERUMPLASMA	FZ
IDENTIFIED DRUGS	840M	RED	SERUM	
IGA SUBCLASS 1	4495L	RED SST	SERUM	FZ
IGA SUBCLASS 2	4496L	SST RED	SERUM	RF
IGF-1	198L	SST	SERUM	RF
IGF-2	704L	RED	SERUM	FZ
IGF-BP1	4296L	RED SST LAV GR	SERUM PLASMA	FZ
IGF-BP3	197L	SST	SERUM	RF
IGG SUBCLASS 1	996L	RED SST	SERUM	RF
IGG SUBCLASS 2	997L	RED SST	SERUM	RF
IGG SUBCLASS 3	998L	RED SST	SERUM	RF
IGG SUBCLASS 4	999L	RED SST	SERUM	RF
IMIPRAMINE (TOFRANIL), SERUM	2001L	RED,LAV,GREEN	SERUM/PLASMA	RF
IMIPRAMINE, SERUM	2002L	RED,LAV,GREEN	SERUM/PLASMA	RF
IMMIGRATION PNL (1SST)	7571	SST	SERUM	
IMMOTILE	131L	SC	SEMEN	RM
IMMUNOFIXATION, 24 HOUR URINE	113L	24 HR URINE	24 HR URINE	RF
IMMUNOFIXATION, SERUM	250L	SST	SERUM	RF
IMMUNOFIXATION, SERUM	250Q	RED SST	SERUM	RF
IMMUNOFIXATION, URINE	899L	SC	URINE	RF
IMMUNOFIXATION, URINE	899Q	SC OR 24HR UC	URINE	RF
IMMUNOGLOBULIN A	628Q			
IMMUNOGLOBULIN A SUBCLASS 1&2	4497L	RED SST	SERUM	RF
IMMUNOGLOBULIN A, SERUM	628L	SST	SERUM	RF
IMMUNOGLOBULIN E, TOTAL	061L	SST	SERUM	RF
IMMUNOGLOBULIN G	630Q			
IMMUNOGLOBULIN G SUBCLASS 1-4	995L	RED SST	SERUM	RF
IMMUNOGLOBULIN G, SERUM	630L	SST	SERUM	RF
IMMUNOGLOBULIN M	629Q			
IMMUNOGLOBULIN M, SERUM	629L	SST	SERUM	RF
IMMUNOGLOBULIN PNL (A,M,G)	4128Q			
IMMUNOGLOBULINS PANEL	4128L	RED SST	SERUM	RF
INFLAMMATORY BOWEL DISEASE	828L	RED SST	SERUM	RF
INFLUENZA A Ag, EIA	3028L	VTM	SWAB	RF
INFLUENZA A AND B CULTURE, RAPID	159Q	SWAB SC VTM	SWAB FLUID	RT
INFLUENZA A H1N1, RT PCR	5097Q	SC SWAB VTM	SWAB FLUID	RF
INFLUENZA A&B Ag, EIA	3030L	VTM	SWAB	RF
INFLUENZA A/B RT PCR, RFLX H1N1	5093Q	SC SWAB VTM	SWAB FLUID	RF
INFLUENZA E A/B Abs QUANT	4376L	RED SST	SERUM	RF
INFLUENZA E B, IgG	4303L	RED SST	SERUM	RF
INITIAL PROFILE	2201	1SST,L	SERUM/PLASMA	RM
INITIAL VISIT NEW PROFILE	7231	2LAV,2SST,1RED	WB,SERUM	RF
INR	910J	BLUE	PLASMA	RF
INR	930B	BLUE	PLASMA	RF
INR	INR	BT	PLASMA	RM
INR	INRL	BLUE	PLASMA	RF
INSULIN	337	SST	SERUM	RF
INSULIN AUTOANTIBODY	4271L	RED SST GR LAV	SERUM PLASMA	FZ
INSULIN DEPENDENT	453L	SST	SERUM	RF
INSURANCE PNL (SST,LAV,UA)	7570			
INTERLEUKIN-6 (IL-6)	3074L	RED SST	SERUM	FZ
INTERPRETATION	0043			
INTERPRETATION	964			
INTERPRETATION	9099			
INTERPRETATION	489L	RED SST		

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
INTERPRETATION	767L	SST	SERUM	RF
INTERPRETATION	86L			
INTRINSIC FACTOR Ab	5211L	RED	SERUM	RF
IRINOTECAN TOXICITY (UGT1A1)	4676L	ACD LAV RT	WHOLE BLOOD	RT
IRON	145	SST	SERUM	RF
IRON AND TIBC	290	SST	SERUM	RF
IRON SATURATION	979	SST	SERUM	RF
IRON SATURATION	%FESA	SST	SERUM	RF
IRON SATURATION	158L	SST	SERUM	RF
IRON, SERUM	145L	SST	SERUM	RF
ISLET CELL DYSFUNCTION I	4262L	RED SST GR LAV	SERUM PLASMA	FZ
ISOPROPANOL, URINE	695L	SC	URINE	RF
JACKSON PANEL-SST,LV	39160	SST,LAV	SERUM, BLOOD	
JAK2 EXON 12 MUTATION ANALYSIS	4600L	LAV GREEN	WHOLE BLOOD	RT
JAK2 V617F MUTATION	3055L	LAV	BLOOD MARROW	RT
Jo-1 Ab	448	SST	SERUM	RF
jose 1	JOSE1			
jose 2	JOSE2			
jose 3	JOSE3			
K033 OAK (WOOD DUST)	4368L	SST	SERUM	RF
K082 LATEX	2038L	SST	SERUM	RF
KEPPRA (LEVETIRACETAM)	470L	RED	SERUM	RF
KETAMINE SCREEN, URINE	871M	SC	URINE	RF
KETONES	U05	SC	URINE	RF
KOH	114L	SWAB	FLUID	RM
LACTIC ACID, PLASMA	514L	GRAY	PLASMA	RF
LACTOFERRIN	4059L	SC	STOOL	FZ
LAKE OB PANEL	2014	2SW,PRO,SST,2L	P-S-FLUID	RM
LAKEPARK FEMALE HORMONE	P120	2SST	SERUM	RM
LAKEPARK GENERAL HEALTH	P122	1SST,L,U	S-P-U	RF
LAKEPARK MALE HORMONE	P121	2SST	SERUM	RF
LAMOTRIGINE, SERUM	205L	RED	SERUM	RT
LAP SCORE	596L	LAV	SLIDES	RM
LD ISOENZYMES	4576L	RED SST	SERUM	RF
LDH	146	SST	SERUM	RF
LDH	4578L	RED SST	SERUM	RF
LDH ISOENZYMES	4576Q	SST RED	SERUM	RF
LDL (calc)	4044L	SST	SERUM	RF
LDL (CALCULATED)	LDL	SST	SERUM	RF
LDL calc	4078L	SST	SERUM	RF
LDL CHOLESTEROL, DIRECT	2031L	SST	SERUM	RF
LDL/HDL RATIO	4045L	SST	SERUM	RF
LDLC	LDLCH	SST	SERUM	RF
LDL-R (REAL)-C	5051L	RED SST	SERUM	RF
LEAD, BLOOD	147L	ROYAL BLUE	WHOLE BLOOD	RF
LEAD, BLOOD	147Q	ROYAL BLUE WB	WHOLE BLOOD	RT
LEAD, BLOOD (PEDIATRIC)	148L	LAV/RB	WHOLE BLOOD	RF
LEAD, URINE	657L	SC	URINE	RF
LEAD/PROTOPORPHYRIN, PEDIATRIC	152L	LAV/RB	WHOLE BLOOD	RT
LEGIONELLA PNEUMONPHILA Abs	4634L	RED SST	SERUM	RF
LEPTIN, SERUM	4038L	RED SST LAV	SERUM PLASMA	RF
LEPTOSPIRA Ab, IgM	262L	SST	SERUM	RF
LEUCINE AMINOPEPTIDASE	342L	RED	SERUM	RF
LEUKEMIA/LYMPHOMA IMMUNOPHENOTYPING	1095L	GR,LAV OR ACD	WB OR BM	RT
LEUKEMIA/LYMPHOMA MONITOR PROFILE	1089L	RED SST	SERUM	RF
LEUKOCYTES STOOL	S12	SC	STOOL	RM
LEUKOCYTES URINE	U12	SC	URINE	RF
LH	102	SST	SERUM	RF
LH	106L	SST	SERUM	RF
LIDOCAINE(XYLOCAINE), SERUM	739L	RED	SERUM	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
LIFE RECOVERY URINE PNL-SC	80010	URINE	URINE	
LIGHTHOUSE COTTAGE UDS	70740	SC	URINE	
LIPASE	148	SST	SERUM	RF
LIPID ASSOC SIALIC ACID	1088L	RED,SST	SERUM	RF
LIPID PANEL	1017	SST	SERUM	RF
LIPID PANEL	1017H	SST	SERUM	RF
LIPID PANEL	4046L	SST	SERUM	RF
LIPID PANEL WITH LDL/HDL RATIO	4047L	SST	SERUM	RF
LIPOPROFILE NMR	352L	LAV,1RED	SERUM	RF
LIPOPROTEIN (a)	620L	SST	SERUM	RF
LIPOPROTEIN A2	4615L	RED SST GR LAV	SERUM PLASMA	FZ
LIPOPROTEIN ELECTROPHORESIS, NMR	352Q	RED SST	SERUM	RF
LIPOPROTEIN PHENOTYPING PANEL	4075L	SST	SERUM	RF
LIPOSUCTION PANEL-2SST,L,B	39240	2SST,LAV,B	BLOOD,PLASM	
LIQUEFACTION TIME	124L	SC	SEMEN	RM
LITHIUM	075	RED	SERUM	RF
LITHIUM	075L	RED	SERUM	RF
LITHIUM	075MM	RED	SERUM	RF
LITHIUM PNL/SST	38342	SST,	SERUM,PLASMA	
LIVER-KIDNEY MICROSOMAL Ab	202L	SST	SERUM	RF
LORAZEPAM (ATIVAN)	796L	RED	SERUM	RF
LORAZEPAM (ATIVAN)	946M	RED GREEN	SERUM PLAMSA	RF
LSD, URINE	88L	SC	URINE	RF
LUPUS ACTIVITY PANEL	4020	2SST,2LAV,1U	S.P.U	RF
LUPUS ANTICOAGULANT	299L	BLUE	CITRATED PLA	FR
LUPUS ANTICOAGULANT	299Q	BLUE FROZEN	PLASMA	FZ
LUPUS ERYTHEMATOSUS	051	RED	SERUM	RF
LUPUS PANEL	4010	3SST,RED,LAV,1U	S-P	RF
LYME AB IgG INTERPRETATION	2016L	RED OR SST	SERUM	RT
LYME Ab, QUANT, IgM	3066L	SST	SERUM	RF
LYME Ab/TOTAL IMMUNOGLOBULINS	632L	SST	SERUM	RF
LYME Abs, TOTAL	636Q	RED SST	SERUM	RF
LYME Abs, WESTERN BLOT	2005Q	RED SST	SERUM	RF
LYME IgG WESTERN BLOT	2005L	RED SST	SERUM	RF
LYME IgG/IgM Ab	636L	RED	SERUM	RF
LYME IGM WB INTERPRETATION	2025L	RED OR SST	SERUM	RT
LYME IgM WESTERN BLOT	2017L	RED OR SST	SERUM	RT
LYMPH ABS	8016B	LAV	WB	RF
LYMPH%	8011B	LAV	WB	RF
LYMPHOCYTE	LYM#J	LAV	WB	RF
LYMPHOCYTE	LYM%J	LAV	WB	RF
LYMPHOCYTE ACTIVATION	1098L	GREEN	WHOLE BLOOD	RT
LYMPHOCYTES	LYMP%	LAV	WB	RF
LYMPHOCYTES #	LYMP	LAV	WB	RF
LYSOZYME, SERUM	2027L	RED OR SST	SERUM	FZ
M ADM-SST,L,SC	70351	SST,LAV,SC	SERUM,PLASM	RM
M. PNEUMONIAE Ab PANEL	397L	SST	SERUM	RF
M. PNEUMONIAE IgG	766L	RED	SERUM	RF
M. PNEUMONIAE IgM	768L	RED	SERUM	RF
M. TUBERCULOSIS CULTURE	941L	SWAB	FLUID	RM
M. TUBERCULOSIS PCR	943L	SWAB	FLUID	RM
M. TUBERCULOSIS, PCR/CULTURE	940L	SWAB	SWAB	RM
M001 PENICILLIUM NOTATUM	2060L	SST	SERUM	RF
M002 CLADOSPORIUM HERBARUM	2061L	SST	SERUM	RF
M003 ASPERGILLUS FUMIGATUS	2062L	SST	SERUM	RF
M004 MUCOR RACEMOSUS	2063L	SST	SERUM	RF
M005 CANDIDA ALBICANS	4228L	SST	SERUM	RF
M006 ALTERNARIA TENUIS	2064L	SST	SERUM	RF
M009 FUSARIUM MONILIFORME	4230L	SST	SERUM	RF
M010 STEMPHYLIUM BOTRYOSUM	2065L	SST	SERUM	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
M012 AUREOBASIDIUM PULLULANS	4231L	SST	SERUM	RF
M013 PHOMA BETAE	4232L	SST	SERUM	RF
M014 EPICOCUM PURPURASCENS	4233L	SST	SERUM	RF
M045 HORMODENDRUM HORDEI	4236L	SST	SERUM	RF
M080 HELMINTHOSPORIUM SPP.	4229L	SST	SERUM	RF
M088 STEMPHYLIUM SOLANI	4234L	SST	SERUM	RF
MACRO	MACR			
MAG Ab, IgM	4535L	SST RED	SERUM	RF
MAGNESIUM	150	SST	SERUM	RF
MAGNESIUM, 24 HR URINE	4503L	24UC ACIDIFIED	URINE (ACID)	RF
MAGNESIUM, RBC	589L	LAV	RBC's	RF
MAIL	MAIL			
MALARIA SMEAR	433L	SLIDE	SMEAR	RM
MALE >40-2SST,L,SC,HCARDS	31033			
Male Admit Panel - 2 sst + lav	30891	2sst lav		
MALE ADMIT PANEL- 2SST + 1LAV	1060	SST,LAV	SERUM,PLASMA	RM
MALE HORMONE-2SST	95403		SERUM	
MALE PANEL- SST+LAV+SC	30911	SST,LAV,SC		
MALE PANEL-LV,B,4SST	71040	LV,B,4SST	ALL	
MANGANESE, BLOOD	167L	ROYAL BLUE	WHOLE BLOOD	RM
MARTIN CTY OXYCODONE	2035	SC	URINE	RF
MARTIN CTY SHERIFF DRUG SCREEN	2030	SC	URINE	RF
MATERNAL AGE AT EDD	450L	SST	SERUM	RF
MATT DRG SCR PNL	39571	SC	URINE	
McCOY PANEL 1	3041		PANEL	ALL
McCOY PANEL 2	3042		PANEL	ALL
McCOY PANEL 3	3043		PANEL	ALL
McCOY PANEL 4	3044		PANEL	ALL
MCH	8005B	LAV		
MCH	MCH	LAV	BLOOD	RM
MCH	MCHJ	LAV	WB	RF
MCHC	8006B	LAV		
MCHC	MCHC	LAV	WHOLE BLOOD	RF
MCHC	MCHCJ	LAV	WB	RF
MCV	8004B	LAV		
MCV	MCV	LAV	PLASMA	RF
MCV	MCVJ	LAV	WB	RF
MDT-3-002CMP,AMY.CA,LIP,GGT,URIC ACID,K	731	LAV	PLASMA	RF
MEDICAL DIRECTOR	0046			
MELATONIN	4625L	FZ LAV PLASMA	PLASMA	FZ
MERCURY, BLOOD	348L	ROYAL BLUE	WHOLE BLOOD	RF
MERCURY, URINE	661L	URINE	URINE	RF
MESORIDAZINE	770L	RED	SERUM	RF
METANEPHRINE, 24 HR URINE	842L	24 HOUR URINE	URINE	RM
METANEPHRINE, PLASMA	4190L	LAV	PLASMA	RF
METANEPHRINE, URINE	841L	SC	URINE	RM
METANEPHRINES FRACT, URINE 24 HOUR	843L	UC	URINE	RF
METANEPHRINES, FRAC., PLASMA, FREE	4187L	LAV	PLASMA	RF
METHADONE	504	SC	URINE	RM
METHADONE	504Q	SC	URINE	RF
METHADONE CONF (GC/MS)	504QU	SC	URINE	RF
METHADONE CONFIRMATION	504M	SC	URINE	RF
METHADONE CONFIRMATION, URINE	4012L	SC	URINE	RF
METHADONE METABOLITE	534Q	UC	URINE	RF
METHAMPHETAMINE	803M	SC	URINE	RF
METHANOL, URINE	693L	SC	URINE	RF
METHAQUALONE	505	SC	URINE	RM
METHAQUALONE	505Q	SC	URINE	RF
METHAQUALONE CONF (GC/MS)	505QU	SC	URINE	RF
METHAQUALONE CONFIRMATION	505M	SC	URINE	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
METHAQUALONE CONFIRMATION, URINE	4013L	SC	URINE	RF
METHODONE, SERUM	912M	RED GREEN	SERUM PLASMA	RF
METHOTREXATE, SERUM	258L	RED	SERUM	RF
METHYLENEDIOXYAMPHETAMINE	804M	SC	URINE	RF
METHYLENEDIOXYMETHAMPHETAMINE	805M	SC	URINE	RF
METHYLMALONIC ACID, SERUM	1161L	RED	SERUM	RF
METHYLMALONIC ACID, URINE	3012L	SC	URINE	RF
METHYLPHENIDATE (RITALIN)	5011L	RED,LAV, GREEN	SERUM,PLASMA	RF
MICRO	MICR			
MICROALBUMIN / CREATININE RATIO	225	URINE	URINE	RM
MICROALBUMIN, 24 HR URINE	842	24 HR URINE	URINE	RF
MICROALBUMIN, URINE	122	SC	URINE	RM
MICROSCOPIC/ADEQUACY EVAL.	0044			
MILES	2			
MILLENNIUM FEMALE - 3 SST 1 LAV	90652			
MILLENNIUM MALE - 3 SST 1 LAV	90651			
MIRTAZAPINE (REMERON), SERUM	860M	RED GREEN	SERUM PLASMA	RF
MIRTAZAPINE (REMERON), URINE	861M	URINE	URINE	RF
MISCELLANEOUS SEND OUT TEST	MISC			
MISSING INFORMATION	4000			
MITOCHONDRIAL (M2) Ab	200Q	RED SST	SERUM	RF
MITOCHONDRIAL(M2) Ab	200L	SST	SERUM	RF
MMR PROFILE	MMR	SST	SERUM	RF
MONO ABS	8017B	LAV	WHOLE BLOOD	RF
MONO QUALITATIVE	052L	SST	WB	RF
MONO TEST	052	SST	SERUM	RF
MONO TEST	052MM	SST RED	SERUM	RF
MONO%	8012B	LAV	WHOLE BLOOD	RF
MONOCYTE	MON#J	LAV	WB	RF
MONOCYTE	MON%J	LAV	WB	RF
MONOCYTES	MONO%	LAV	WHOLE BLOOD	RF
MONOCYTES #	MONO	LAV	WHOLE BLOOD	RF
MORPHINE, URINE QUANT	872M	SC	URINE	RF
MOTILITY	127L	SC	SEMEN	RM
MPV	8008B	LAV	WHOLE BLOOD	RF
MPV	MPV	LAV	WHOLE BLOOD	RF
MTHFR, DNA ANALYSIS	4127L	SWAB LAV ACD	WB SWAB	RT RF
MUCIN CLOT TEST	351L	GN,LAV	SYNOVIAL FLU	RF
MUCUS THREADS	U21	SC	URINE	RF
MULTI-CARE PANEL	30771	2SST,LAV, SC	SRM,PLAS,URI	
MULTIPLE GESTATION	454L	SST	SERUM	RF
MUMPS Ab, IgG	396	SST	SERUM	RF
MUMPS Abs, IgM	395L	SST	SERUM	RF
MURTHY ANEMIA	36021	SST,LAV	SERUM/WHOLE	RF
MURTHY ARTHRITIS	36022	SST,LAV	SERUM/WHOLE	RF
MUSK ANTIBODY TITER	4225L	RED	SERUM	RF
MYASTHENIA GRAVIS PNL III	4324L	RED SST	SERUM	RF
MYCOPHENOLIC ACID & METAB	4696L	RED LAV	SERUM/PLASMA	RF
MYCOPLASMA/UREAPLASMA	984L			
MYCOPLASMA/UREAPLASMA CULTURE	984Q	SC VTM	FLUID	RF
MYELO	MYEL	LAV	WHOLE BLOOD	RF
MYOGLOBIN, SERUM	353L	SST	SERUM	RF
MYOGLOBIN,URINE	4178L	UC	URINE	RF
MYOSITIS PANEL	4635L	RED SST	SERUM	RF
N. MENINGITIDIS SEROGROUP PANEL IgG	4048L	SST	SERUM	RF
NAGEL PROFILE 1	1160	SST,L,S.C.	S-P-U	RM
NAGEL PROFILE 2	1160A	SST,L,S.C.,STOO	S-P-U-STOOL	RM
NAGEL PROFILE 3	1160B	SST,L,S.C.,STOO	S-P-U-STOOL	RM
NAPA, SERUM	023L	RED	SERUM	RF
NASAL CULTURE	611	SWAB	SWAB	RM

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
NATURAL KILLER CELL PNL	4307L	LAV ACD	WHOLE BLOOD	RT
NEISSERIA GONORRHOEAE	416	PROBETEC	FLUID	RM
NEISSERIA GONORRHOEAE, DNA	612Q	THIN PREP	THIN PREP	RT
NEISSERIA GONORRHOEAE, NAA	4400L	URINE OR SWAB	URINE/SWAB	RF
NEUROFIBROMATOSIS TYPE 1	4191L	GREEN	WHOLE BLOOD	RT
NEURON SPECIFIC ENOLASE	2087L	RED OR SST	SERUM	RF
NEUROPATHY PROFILE II	4624L	2SST	SERUM	RF
NEUT ABS	8015B	LAV	WB	RF
NEUT%	8010B	LAV	WHOLE BLOOD	RF
NEUTROPHIL	NEU#J	LAV	WB	RF
NEUTROPHILS	NEU%J	LAV	WB	RF
NEUTROPHILS	SEGM%	LAV	WHOLE BLOOD	RF
NEUTROPHILS #	SEGM	LAV	WHOLE BLOOD	RF
NEW PATIENT I	4030	SST,LAV	S,P	RF
NEW PATIENT II	4040	3 SST	S	RF
NEW PATIENT III	4050	2 SST	S	RF
NEW PATIENT IV	4060	SST	S	RF
NEXSLIM PANEL	31790	SST	SERUM	RF
NEXT VISIT NEW PROFILE	7232	GRAY,RED	SERUM	RF
NHTC DEPAKOTE PNL-LAV*FRZ LAV*SST	38431	2LAV, SST	SERUM,PLASMA	
NICKEL, PLASMA	168L	ROYAL BLUE	WHOLE BLOOD	RM
NICOTINE QUANTITATION	846M	SC	URINE	RF
NICOTINE/COTININE CONFIRMATION	844M		URINE	RF
NICOTINE/COTININE SCREEN	845M	SC	URINE	RF
NICOTINIC ACID (B3)	4220L	LAV	PLASMA(LIGHT	FZ
NIEMANN-PICK, TYPE A&B, DNA ANALYSIS	4550L	LAV ACD	WHOLE BLOOD	RF RM
NITRITE	U11	SC	URINE	RM
non HDL CHOL	5000H	SST	SERUM	RF
NORBUPRENORPHINE SCREEN	854M	SC	URINE	RF
NOREPINEPHRINE, 24HR URINE	137L	24 HR URINE	URINE	RF
NOREPINEPHRINE, URINE	136L	SC	URINE	RF
NORMAL MORPHOLOGY	135L	SC	SEMEN	RM
NORMETANEPHRINE, 24 HR URINE	350L	24 HR URINE	URINE	RF
NORMETANEPHRINE, PLASMA	4189L	LAV	PLASMA	RF
NORMETANEPHRINE, URINE	349L	SC	URINE	RF
NOROVIRUS, RT-PCR	696L	SC	STOOL	RF
NORTRIPTYLINE, SERUM	866L	RED/LAV/GR	SERUM/PLASMA	RF
NRBC	NRBC	LAV	WHOLE BLOOD	RF
N-TELOPEPTIDE, SERUM	588L	SST	SERUM	RF
N-TELOPEPTIDE, URINE	499L	SC	URINE	RF
O201 TOBACCO DUST	4237L	SST	SERUM	RF
OB PANEL	2015	1SST,R,2LAV	SERUM,PLASMA	RM
OB PANEL TSINKER & LOSS	912	1SST,R,2LAV	SERUM,PLASMA	RM
OBSTETRIC PANEL 1	4043	1SST,R,2LAV	SERUM,PLASMA	RM
OBSTETRIC PANEL, DUMENY	3011	1SST,R,2LAV	SERUM,PLASMA	RM
OBSTETRICS PANEL	1018	1SST,R,2LAV	SERUM,PLASMA	RM
OCCULT BLOOD	031	SC	STOOL	RM
OLANZAPINE (ZYPREXA)	976M	RED GREEN	SERUM PLASMA	RF
OPIATE ANALYSIS, SERUM	960M	RED SST	SERUM	RF
OPIATE CONF (GC/MS)	509QU	SC	URINE	RF
OPIATE CONF (GC/MS) EXPANDED	969Q	SC	URINE	RF
OPIATE CONF. (MEDTOX 9376)	969M	SC	URINE	RF
OPIATE CONFIRMATION	509M	SC	URINE	RF
OPIATE CONFIRMATION, BLOOD	3097L	GRAY	WHOLE BLOOD	RF
OPIATE CONFIRMATION, URINE	4005L	SC	URINE	RF
OPIATES	509	SC	URINE	RF
OPIATES	509Q	SC	URINE	RF
ORGANIC ACID INTERPRETATION	5027L	TT	URINE	FZ
ORGANIC ACID, PLASMA	3001L	GREEN	PLASMA	FZ
ORGANIC ACID, URINE PANEL	5026L	TT	URINE	FZ

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
ORGANISM IDENTIFICATION	625			
OSBR RISK	462L	SST	SERUM	RF
OSMOLALITY, SERUM	358L	SST	SERUM	RF
OSMOLALITY, URINE	359L	SC	URINE	RF
OSMOLALITY, URINE 24 HOUR	360L	24 HR URINE	URINE	RF
OSTEOCALCIN, SERUM	5033L	RED SST	SERUM	FZ
OTHER MICROSCOPIC OBSERVATIONS	134L	SC	SEMEN	RM
OVA AND PARASITE	032L	STOOL	STOOL	RM
OVA AND PARASITES	032Q	OVP KIT	STOOL	RF
OVARIAN PANEL	2020	2SST	SERUM	RF
OXALATE, 24 HR URINE	4504L	24UC ACIDIFIED	URINE (ACID)	RF
OXAZEPAM, QUANT GC/MS	832M	UC	URINE	RF
OXCARBAZEPINE (TRILEPTOL)	305L	RED	SERUM	RF
OXYCODONE	171	SC	URINE	RF
OXYCODONE	171Q	SC	URINE	RF
OXYCODONE	171QU	SC	URINE	RF
OXYCODONE	529M	SC	URINE	RF
OXYCODONE CONFIRMATION	510M	SC	URINE	RF
OXYCODONE CONFIRMATION, URINE	4014L	SC	URINE	RF
OXYCODONE, BLOOD	968M	RED GREEN	SERUM PLASMA	RF
OXYCODONE, QUANT	838M	RED	SERUM	RF
P PARTNERS ADM-SST,LV	70901	LV,SST	SERUM,BLOOD	
P PARTNERS UDS	70630	SC	URINE	
p18 AB	2015L	RED OR SST	SERUM	RT
p18 Ab	5232L	SST	SERUM	RF
p23 AB	2014L	RED OR SST	SERUM	RT
p23 AB	2020L	RED OR SST	SERUM	RT
p24 Ab	5231L	SST	SERUM	RF
p28 AB	2013L	RED OR SST	SERUM	RT
p30 AB	2012L	RED OR SST	SERUM	RT
p31 Ab	5230L	SST	SERUM	RF
p39 AB	2011L	RED OR SST	SERUM	RT
p39 AB	2019L	RED OR SST	SERUM	RT
p41 AB	2010L	RED OR SST	SERUM	RT
p41 AB	2018L	RED OR SST	SERUM	RT
p45 AB	2009L	RED OR SST	SERUM	RT
p51 Ab	5227L	SST	SERUM	RF
p55 Ab	5233L	SST	SERUM	RF
p58 AB	2008L	RED OR SST	SERUM	RT
p65 AB	5226L	SST	SERUM	RF
p66 AB	2007L	RED OR SST	SERUM	RT
p93 AB	2006L	RED SST	SERUM	RT
PAIN MGMT PANEL- URINE	39001		URINE	RF
PALMS WELLINGTON PRE-OP-SST,L,B	71080	SST,LAV,B	BLD,URIN,PL	
PANCREATIC ELASTASE, FECAL	1045L	SC	STOOL	FZ
PANCREATIC POLYPEPTIDE	5234L	RED SST	SERUM	FZ
panel 100	1111			
panel firm	2222			
PAP SMEAR CONVENTIONAL	004	SLIDE	SLIDE	RM
PAP SMEAR, CONVENTIONAL	004Q	SLIDE	PAP	RT
PARANEOPlastic PROFILE 1	4567L	RED SST	SERUM	RF
PARASITE EXAM, BLOOD	486L	LAV SMEAR	WB SMEAR	RF
PARASITE ID, WORM	4715L	SC		
PARATHYROID HORMONE-INTACT	362	RED SST	SERUM	RF
PARIETAL CELL ANTIBODY	642L	SST	SERUM	RF
PAROXETINE (PAXIL)	673	RED	SERUM	RF
PARVOVIRUS B19, IgG	478L	SST	SERUM	RF
PARVOVIRUS B19, IgM	479L	SST	SERUM	RF
PATHOLOGY REVIEW	5555	SLIDE	SLIDE	RM
PBI COMP DRUG PROFILE- SC	71100	URINE	URINE	

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
PERCENT FREE CORTISOL	700L	SST	SERUM	RF
PERIPHERAL SMEAR REVIEW	7411	SLIDE	SLIDE	RM
PERPHENAZINE (TRILAFON)	279M	RED SST	SERUM	RF
PESTICIDE EXPOSURE PANEL	4701L	RED GREEN	SERUM PLASMA	RT
pH	535	UC	URINE	RF
pH	126L	SC	SEMEN	RM
pH	U08	SC	URINE	RF
PHA-CD4 ACTIVATION	1097L	GREEN	WHOLE BLOOD	RT
PHENCYCLIDINE	508	SC	URINE	RF
PHENCYCLIDINE	508Q	SC	URINE	RF
PHENCYCLIDINE CONF (GC/MS)	508QU	SC	URINE	RF
PHENCYCLIDINE CONFIRMATION	508M	SC	URINE	RF
PHENOBARBITAL	076	RED	SERUM	RF
PHENOBARBITAL	076L	RED	SERUM	RF
PHENOBARBITAL	076MM	RED	SERUM	RF
PHENYTOIN	074MM	RED	SERUM	RF
PHENYTOIN (DILANTIN)	074	RED	SERUM	RM
PHENYTOIN (DILANTIN)	074H	RED	SERUM	RF
PHENYTOIN, FREE	020L	RED	SERUM	RF
PHENYTOIN, FREE AND TOTAL	4515L	RED SST	SERUM	RF
PHENYTOIN, TOTAL	4516L	RED SST	SERUM	RF
PHLEBOTOMY CHARGE	0007			
PHOENIX DRUG PANEL 10	81710		URINE	
PHOSPHOLIPIDS, SERUM	654L	RED	SERUM	RF
PHOSPHORUS	151	SST	SERUM	RF
PHOSPHORUS, 24HR URINE	1094L	24 HR UC 6NHCL	URINE	RT
PHOSPHORUS, URINE	1093L	SC	URINE	RT
PHY WT.LOSS PNL-SST,LAV	31361	SST,LAV	BLOOD,SERUM	
PHYSICAL PNL- 2SST,LAV,SC	70780	SST,LAV,SC	BLD,PLS,UR	
PINWORM IDENTIFICATION	3015Q	SC TAPE PREP	TAPE PREP	RT
PINWORM PREP, ENTEROBIUS	3015L	GLASS SLIDE	TAPE TEST	RT
PLAC TEST	4224L	SST	SERUM	RF
PLASMINOGEN ACTIVITY	5037L	BLUE	PLASMA	FZ
PLASMINOGEN ANTIGEN	4116L	BLUE	PLASMA	FZ
PLASMODIUM sp. DNA by PCR	886L	LAV or ACD	WHOLE BLOOD	RT or RF
PLATELET	033	LAV	WHOLE BLOOD	RF
PLATELET	8007B	LAV	WHOLE BLOOD	RF
PLATELET ANTIBODY, DIRECT	4280L	ACD LAV	WHOLE BLOOD	RT
PLATELET ANTIBODY, SERUM	692L	SST	SERUM	RF
PLATELETS	033J	LAV	WB	RF
PLATELETS	033MM	LAV	WHOLE BLOOD	RF
PLT ESTIMATE	PLTE	LAV	WHOLE BLOOD	RF
PLYCHROM	PLYC	LAV	WHOLE BLOOD	RF
PNEUMOCOCCAL IMMUNITY (14)	773L	SST	SERUM	RF
PNH EVALUATION	4277L	ACD GREEN	WHOLE BLOOD	RT
POIK	POIK	LAV	WHOLE BLOOD	RF
POLIOVIRUS Ab	3031L	RED SST	SERUM	RF
PORPHOBILINOGEN, URINE	584L	SC	URINE	RF PROTECT FROM LIGHT
PORPHYRINS QUANT, 24 HOUR URINE	720L	24 HR URINE	URINE 24 HR	RF
PORPHYRINS QUANT, URINE RANDOM	735L	SC	URINE	RF
POTASSIUM	152	SST	SERUM	RF
POTASSIUM	1521	SST	SERUM	RF
POTASSIUM	152H	SST	SERUM	RT
POTASSIUM	152J	SST	SERUM	RF
POTASSIUM	6001B	SST	SERUM	RF
POTASSIUM, 24 HR URINE	516L	24 HR URINE	URINE	RF
POTASSIUM, URINE	515L	SC	URINE	RF
POWERS PROFILE 1 (SST,LAV,BL)	90811			
POWERS PROFILE 2 (2SST,LAV,BL)	90812			
POWERS PROFILE 3 (2SST,LAV,BL)	90813			

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
POWERS PROFILE 4 (2SST,2LAV,BL)	90814			
POWERS THYROID PROFILE (2SST)	90815			
PREALBUMIN	369MM	RED SST	SERUM	RF
PRE-ALBUMIN	369	SST	SERUM	RF
PRE-BETA	4084L	SST	SERUM	RF
PREGABALIN (LYRICA)	859M	RED	SERUM	RF
PREGABALIN, SERUM	5031L	RED	SERUM	RF
PREGNENOLONE	480L	RED SST LAV FZ	SERUM PLASMA	FZ
PRE-NATAL PROFILE I	1067	2L,R,SST,SC,PRO	P-S-U-F	RM
PRE-NATAL PROFILE II	1068	SST,R,L,S.C.	P-S-U	RM
PRE-OP PANEL	3006	L,B,SC,SST	P-S-U	RM
PRE-OPERATIVE PANEL	1004	L,B,SST,S.C.	S-P-U	RM
PRILOCAINE	740L		PLASMA	RF
PRIMIDONE	651L	RED	SERUM	RF
PRO	PROT	SC	URINE	RF
PROCAINAMIDE, SERUM	022L	RED	SERUM	RF
PROCALCITONIN	4559L	RED SST	serum	RF
PROFILE 1 (2SST, 1LAV, 1UA)	70161			
PROFILE 2 (2SST, 1LAV, 1UA)	70162			
PROGESTERONE	162	SST	SERUM	RM
PROGESTERONE	162L	SST	SERUM	RF
PROGRAF (TACROLIMUS)	1026L	LAV,GREEN	PLASMA	RF
PROINSULIN	314L	RED	SERUM	FZ
PROINSULIN AS PERCENT	4264L	SST RED GR LAV	SERUM PLASMA	FZ
PROLACTIN	103	SST	SERUM	RF
PROLACTIN	107L	SST	SERUM	RF
PROPAFENONE, SERUM	953L	RED	SERUM	RF
PROPOXYPHENE	507	SC	URINE	RF
PROPOXYPHENE	507Q	SC	URINE	RF
PROPOXYPHENE CONF (GC/MS)	507QU	SC	URINE	RF
PROPOXYPHENE CONFIRMATION	507M	SC	URINE	RF
PROPOXYPHENE, BLOOD	967M	RED GREEN	SERUM PLASMA	RF
PROSTAGLANDIN e2, URINE	4645L	URINE	URINE	FZ
PROSTATIC ACID PHOS	913L	RED OR SST	SERUM	RF
PROTEIN	U09	SC	URINE	RF
PROTEIN C ANTIGEN	857L	BLUE	PLASMA	FZ
PROTEIN C-FUNCTIONAL	294L	BLUE	SERUM	FZ
PROTEIN ELECTROPHORESIS, 24HR URINE	116Q	24HR UC	URINE	RF
PROTEIN ELECTROPHORESIS, URINE RANDOM	112L	SC	URINE	RF
PROTEIN ELECTROPHORESIS+IFE 24HR URINE	113Q	24HR UC	URINE	RF
PROTEIN S Ag, FREE	859L	BLUE	PLASMA	FZ
PROTEIN S Ag, TOTAL	858L	BLUE	PLASMA	FZ
PROTEIN S FUNCTIONAL	863L	BLUE	PLASMA	FZ
PROTEIN, BODY FLUID	4508L	SC	BODY FLUID	RT
PROTHROMBIN IgG Abs	871L	RED OR SST	SERUM	RF
PROTHROMBIN IgM Abs	872L	RED OR SST	SERUM	RF
PROTIME	034B	BLUE	PLASMA	RF
PROTIME	034J	BLUE	PLASMA	RF
PROTIME	900J	BLUE	PLASMA	RF
PROTIME	920B	BLUE	PLASMA	RF
PROTIME	944L	BLUE	PLASMA	RF
PROTOPORPHYRIN FREE (FEP)	149L	LAV/RB	WHOLE BLOOD	RM
PSA PRE-OP PNL(BLU,LAV,UA)	9370	BLU LAV UA	PANEL	RF
PSA TOTAL & FREE	924	SST	SERUM	RF
PSA, DIAGNOSTIC	90841	RED SST	SERUM	RF
PSA, FREE	552	SST	SERUM	RF
PSA, MONITOR	2024	SST	SERUM	RF
PSA, SCREEN	044	SST	SERUM	RF
PSA, TOTAL	044L	SST	SERUM	RF
PSEUDOEPHEDRINE/EPHEDRINE PANEL	892M	SC	URINE	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
PT	018	BLUE	PLASMA	RF
PT	034	BLUE	PLASMA	RF
PT	034H	BLUE	PLASMA	RF
PT	034L	BLUE	PLASMA	RF
PT	034MM	BLUE	PLASMA	RF
PT MIXING STUDY	946L	BLUE	PLASMA	FZ
PTH-rP (RELATED PROTEIN)	4692L	LAV/TYRASYLOL	PLASMA	FZ
PTT	035	BLUE	PLASMA	RF
PTT	035H	BLUE	PLASMA	RF
PTT	035L	BLUE	PLASMA	RF
PTT	035MM	BLUE	PLASMA	RF
PTT	751J	BLUE	PLASMA	RF
PUC PANEL- SST,LV,B	30060	sst, lv , b	serum,plasma	
PYRIDINIUM CROSSLINKS	4399L	SC	URINE	FZ
PYRUVIC ACID, BLOOD	2096L	GRAY	PA SUPER	RF
QUALITATIVE ANALYSIS	HC36			
QUALITY OF MOTION	128L	SC	SEMEN	RM
QUANTATIVE ANALYSIS (ALC TEST PNL)	1073	SC	URINE	RF
QUANTIFERON (QFT)TB	4584L	LABCORP KIT	WHOLE BLOOD	RT
QUANTITATIVE CONFIRMATION	1040			
QUEST SEND OUT	593			
QUETIAPINE (SEROQUEL)	907L	RED/SST	SERUM	RF
QUETIAPINE (SEROQUEL),SERUM	971M	RED GREEN	SERUM PLASMA	RF
QUINIDINE, SERUM	024L	RED	SERUM	RF
RA LATEX TURBID	053L	SST	SERUM	RF
RACE	493	SST	SERUM	RF
RACE	451L	SST	SERUM	RF
RAPID PLASMA REAGIN QUANT	041L	SST	SERUM	RF
RAPIDLY PROGRESSIVE	129L	SC	SEMEN	RM
RBC	8001B	LAV	WHOLE BLOOD	RF
RBC	RBC	LAV	WHOLE BLOOD	RF
RBC	RBCJ	LAV	WB	RF
RBC'S	U14	SC	URINE	RM
RDW	8009B	LAV	WHOLE BLOOD	RF
RDW	RDW	LAV	WHOLE BLOOD	RF
RDW	RDWJ	LAV	WB	RF
RECOVERY DRUG SCREEN	38401	SC	URINE	
RECOVERY FEMALE-2SST-LAV	30922			
RECOVERY MALE-2SST,L	30921	SST/LAV	BLOOD,PLASMA	
RECOVERY PLC UDS	30920		URINE	
RENAISSANCE UDR10	70111	SC	URINE	
RENAL FUNCTION PANEL	1016	SST	SERUM	RF
RENIN, PLASMA	372L	LAV OR ACD	PLASMA	FZ
REPTILASE TIME	4211L	BLUE	PLASMA	FZ
RESPIRATORY INFECT PROF B	2095L	SST RED	SERUM	RT
RETIC COUNT	036	LAV	PLASMA	RF
RETICULIN IgA Ab	647L	SST	SERUM	RF
RETICULIN IgG Ab	648L	SST	SERUM	RF
RETICULOCYTE COUNT	036H	LAV	WHOLE BLOOD	RF
RETICULOCYTE COUNT	036L	LAV	PLASMA	RF
RETICULOCYTE COUNT	036MM	LAV	WHOLE BLOOD	RF
RETT SYNDROME, DNA	1091L	LAV OR YELLOW	WHOLE BLOOD	RT
RH TYPE	404L	LAV	PLASMA	RM
RHEUMATOID FACTOR	053	SST	SERUM	RF
RHEUMATOID FACTOR(QUANT)	439L	SST	SERUM	RM
RHEUMATOID PROFILE	2013	1SST,R,LAV	SERUM,PLASMA	RM
RIBA (HCV CONFIRMATION)	2100L	RED SST	SERUM	RF
RIGHT PLC ADMIT-SST,LAV	31351			
RISPERIDONE AND METABOLITE	953M	RED GREEN	SERUM PLASMA	RF
RISPERIDONE(RISPERDAL),SERUM	4241L	RED	SERUM	RT

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
RITALIN	913M	SC	URINE	RF
RMSF, IgG, EIA	746L	SST	SERUM	RF
RNP ANTIBODIES	447L	SST	SERUM	RF
ROTAVIRUS ANTIGEN DETECTION	1082Q	SC	STOOL	RF
ROTOVIRUS Ag	1082L	SC	STOOL	RF
ROUTINE PANEL (2SST,1LAV)	91511			
RPR	054	SST	SERUM	RF
RSV Ab, QUANT	496L	SST	SERUM	RF
RSV Ag, EIA	497L	VTM	SWAB	RF
RUBELLA Ab, IgG	921	SST	SERUM	RF
RUBELLA IgG	105B	SST	SERUM	RF
RUBELLA, IgM	414L	RED SST	SERUM	RF
RUBEOLA Ab, IgM	410L	SST	SERUM	RF
RUBEOLA ANTIBODY (MEASLES), IgG	398	SST	SERUM	RF
S FLA PAIN PNL-U	70820	SC	URINE	
S.M.CARDIO BLD,fst-2SST,LAV	31430	SST,LAV	SERUM,PLASMA	
S.M.CARDIO HMG-SST	31431	SST	SERUM	
SACCHAROMYCES CEREVISIAE PANEL	302L	SST	SERUM	RF
SALICYLATE, SERUM	025L	RED	SERUM	RF
SALIVARY CORTISOL	4380L	SCT	SALIVA	RF
SALZMAN DRUG PANEL	39321	SC	URINE	RF
SCHISTOSOMA IgG	4330L	RED SST	SERUM TT	RF
SCHISTOSOMIASIS IgG	485L		SERUM	RF
SCLERODERMA-70 Ab	733	SST	SERUM	RF
SCLERODERMA-70 Ab	733L	SST	SERUM	RM
SEDIMENTATION RATE	037	LAV	WHOLE BLOOD	RF
SEDIMENTATION RATE	037B	LAV	WHOLE BLOOD	RF
SEDIMENTATION RATE	037MM	LAV	WHOLE BLOOD	RF
SEIKEL DRUG SCREEN	80090	SC	URINE	RF
SELENIUM, BLOOD	4120L	RB LAV	WHOLE BLOOD	RT
SEMEN ANALYSIS, BASIC	120L	SC	SEMEN	RM
SEMEN ANALYSIS, COMPLETE	120Q	SC	FLUID	RT
SEMEN ANALYSIS, POST VASECTOMY	4591L	SC	SEMEN	RT
SEMI QUANT ANALYSIS	HC37			
SEMI QUANT ANALYSIS (TEST)	8888			
SEMI-QUANT PANEL (TEST)	88888			
SENSITIVITY	2025			
SENSITIVITY STUDY	083			
SENSITIVITY STUDY	244			
SERFER PATRON PANEL I	7000	1SST,1LAV,1U,1B	S,P,U	RF
SEROTONIN	374L	RED	SERUM	RF
SERTRALINE	833M	RED	SERUM	RF
SERUM PROTEIN ELECTROPHORESIS	111L	SST	SERUM	RF
SERUM PROTEIN ELECTROPHORESIS	111Q	RED SST	SERUM	RF
SEX	492	SST	SERUM	RF
SEX HORMONE BINDING GLOB, SERUM	400L	RED	SERUM	RF
SHAWN PROF W/O THYROID	1064	1SST,LAV	SERUM,PLASMA	RM
SHAWN PROF W/THYROID	1063	1SST,LAV,R	SERUM,PLASMA	RM
SHEIK BASIC-2SST,L,SC	31301	SST,LAV,SC	BL,PLAS, UR	
SHEIK FEMALE-SST	31035			
SHIPE PANEL-2LV,2SST,SC	38321	2SST,2LAV,SC	UR,BLD,SERUM	
SICKLE CELL	SICK	LAV	WB	RF
SICKLE CELL SCREEN	039	LAV	WB	RF
SILLER PANEL 5	2042	SST,L,SC	S,P,U	RM
SIROLIMUS (RAPAMUNE)	820L	LAV	WHOLE BLOOD	RF
SKIN CULTURE	657	SWAB	FLUID	RM
SLUGGISH	130L	SC	SEMEN	RM
Sm/RNP Ab	447	SST	SERUM	RF
SMITH (Sm) Ab	570	SST	SERUM	RF
SMOOTH MUSCLE Ab (ACTIN)	755L	SST	SERUM	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
SMOOTH MUSCLE Ab (ACTIN)	755Q	RED SST	SERUM	RF
SODIUM	155	SST	SERUM	RF
SODIUM	155H	SST	SERUM	RF
SODIUM	155J	SST	SERUM	RF
SODIUM	6000B	SST	SERUM	RF
SODIUM, 24 HR URINE	675L	24 HOUR UC	URINE	RF
SODIUM, URINE	375L	SC	URINE	RF
SOLUABLE TRANSFERRIN RECEPTOR	4215L	RED GR LAV	SERUM PLASMA	RF
SPECIFIC GRAVITY	536	UC	URINE	RF
SPECIFIC GRAVITY	U06	SC	URINE	RF
SPECIMEN IDENTITY	4001			
SPHERO	SPHE	LAV	WB	RM
SPINAL MUSCULAR ATROPHY	4678L	LAV ACD STAT	WHOLE BLOOD	RT
SPUTUM CULTURE	610	SC,SWAB	SPUTUM	RM
SPUTUM CYTOLOGY	4320L	SC	SPUTUM	RT
SRP	4245L	SST LAV	SERUM,PLASMA	RF
SS-A/Ro Ab	313	SST	SERUM	RF
SSB/La Ab	711	SST	SERUM	RF
ST LUKES DRUG SCREEN PANEL (8)	3054	UC	URINE	RF
STAT TEST	911			
STONE ANALYSIS	905L	SC	STONE	RF
STONE ANALYSIS	905Q	SC	STONE	RT
STOOL CULTURE	088	SC	STOOL	RM
STOOL CULTURE, YERSINIA	4594L	PARAPAK ORANGE	STOOL	RF
STREP A SCREEN	089	SWAB	FLUID	RM
STREP B SCREEN	1077	SWAB	SWAB	RT
STREP GROUP B SCREEN	1066	SWAB	FLUID	RM
STREP PNEUMO Ag	909L	RED SST	SERUM	RF
STREPTOMYCIN	897L	RED	SERUM	RF
STRONGYLOIDES Ab, IgG	3075L	RED SST	SERUM	RF
SULFATE, 24HR URINE	4124L	24HR UC	24HR URINE	RF
SULFATE, URINE	4125L	UC	URINE	RF
SUNLIGHT ADM PNL- SST,LV	70840	LAV SST	BLOOD/PLASMA	
SURVEILLANCE PROFILE	2205	1SST 2L 1Y	S-P	RM FZ
SYPHILIS IgG	227	SST	serum	rf
SYSTEMIC LUPUS PROFILE A	550L	SST	SERUM	RF
T CELL ACTIVATION, CD8 SUBSETS	3062L	LAV OR ACD	WHOLE BLOOD	RT
T. PALLIDUM Ab (TP-PA)	4388L	RED SST	SERUM	RF
T001 MAPLE/BOX ELDER	4355L	SST	SERUM	RF
T005 BEECH, AMERICAN	4357L	SST	SERUM	RF
T006 CEDAR, MOUNTAIN	3018L	SST	SERUM	RF
T007 OAK, WHITE	3019L	SST	SERUM	RF
T008 ELM, AMERICAN (WHITE)	2066L	SST	SERUM	RF
T010 WALNUT, POLLEN	4358L	SST	SERUM	RF
T016 PINE, WHITE	4419L	SST	SERUM	RF
T018 EUCALYPTUS	2067L	SST	SERUM	RF
T021 MELALEUCA	4420L	SST	SERUM	RF
T027 MAPLE, RED	2068L	SST	SERUM	RF
T030 BIRCH, WHITE	3017L	SST	SERUM	RF
T032 WILLOW, BLACK	4360L	SST	SERUM	RF
T037 CYPRESS, WHITE BALD	2069L	SST	SERUM	RF
T041 HICKORY, WHITE	4361L	SST	SERUM	RF
T061 SYCAMORE, AMERICAN	4359L	SST	SERUM	RF
T070 MULBERRY, WHITE	4325L	RED SST	SERUM	RF
T071 MULBERRY, RED	2081L	SST	SERUM	RF
T072 PALM, QUEEN	2074L	SST	SERUM	RF
T073 PINE, AUSTRALIAN	4421L	SST	SERUM	RF
T083 MANGO TREE	2070L	SST	SERUM	RF
T096 WHITE POPLAR	4362L	SST	SERUM	RF
T103 OAK, LIVE / VIRGINIA	2071L	SST	SERUM	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
T18 (BY AGE)	466L	SST	SERUM	RF
T18 RISK	465L	SST	SERUM	RF
T201 SPRUCE	3020L	SST	SERUM	RF
T202 ALDER, SMOOTH	4356L	SST	SERUM	RF
T210 PRIVET, COMMON	2072L	SST	SERUM	RF
T218 BAYBERRY / SWEET GALE	2073L	SST	SERUM	RF
T3 TOTAL	108	SST	SERUM	RF
T3 UPTAKE	109	SST	SERUM	RF
T3, FREE	440	SST	SERUM	RF
T3, REVERSE	742L	SST	SERUM	RF
T4, FREE	106	SST	SERUM	RF
T4, TOTAL	107	SST	SERUM	RF
TARGET	TARG	LAV	WHOLE BLOOD	RF
TAY SACHS DISEASE LEUKOCYTES	4521L	ACD	WHOLE BLOOD	RF
TAY-SACHS, DNA ANALYSIS	345L	SST	SERUM	RF
TCFW HORMONE PNL-2SST	90301	SST	SERUM	
TEGRETOL PNL/SST,LV	38343	SST,LAV	PLASMA,BLOOD	
TEICHOIC ACID Ab	445L	SST	SERUM	RF
TEMAZEPAM (RESTORIL)	945M	RED GREEN	SERUM PLASMA	RF
TEST CHALLENGES PANEL	1037T			
TEST CODE	TEST			
TEST ENABLED	TEST3			
TEST IN QUESTION	9999			
TEST IN QUESTION	99999			
TEST KITS	0008			
TEST PANEL	TEST1			
TEST PENDING	999			
TESTOSTERONE FREE&WEAKLY BOUND	4536L	RED SST	SERUM	RF
TESTOSTERONE, FREE	377L	SST	SERUM	RF
TESTOSTERONE, TOTAL	110	SST	SERUM	RF
TESTOSTERONE, TOTAL	110L	SST	SERUM	RF
TESTOSTERONE, TOTAL & FREE	377Q	RED SST	SERUM	RF
TETANUS ANTITOXOID IgG	173L	RED OR SST	SERUM	RF
THEOPHYLLINE	077	RED	SERUM	RF
THIN PREP	1027	ThP	FLUID	RM
THIN PREP	1027Q	THIN PREP	PAP	RT
THIN PREP + HPV HIGH	0050Q	THIN PREP	PAP VIAL	RT
THIN PREP + HPV HIGH/LOW	0041Q	THIN PREP	PAP VIAL	RT
THIN PREP RFLX HPV (ASCU OR >)	0042	THIN PREP	PAP	RM
THIN PREP RFLX HPV HIGH	0040Q	THIN PREP	PAP	RT
THIN PREP RFLX HPV HIGH/LOW	0042Q	THIN PREP	PAP	RT
THIOPURINE METHYLTRANSFERASE	4100L	GREEN LAV	WHOLE BLOOD	RF
THIORIDAZINE	769L	RED	SERUM	RF
THROAT CULTURE	087	SWAB	FLUID	RM
THROMBIN TIME	974L	BLUE	PLASMA	RF
THROMBOTIC RISK CONGENITAL	HC1			
THROMBOTIC RISK, CONGENITAL	5039L	BLUE	PLASMA	FZ
THYROGLOBULIN	379	RED SST	SERUM	RF
THYROGLOBULIN Ab	308	RED SST	SERUM	RF
THYROGLOBULIN PANEL	318	RED SST	SERUM	RF
THYROID ANTIBODIES	309	SST RED	SERUM	RF
THYROID FOLLOW UP-SST	95401	SST	SERUM	
THYROID FUNCTION PANEL	1022	SST	SERUM	RF
THYROID NEW PATIENT-2SST	95400	SST	SERUM	
THYROID PANEL (1SST)	36023	SST	SERUM	RF
THYROID PANEL, COMPREHENSIVE	1023	SST	SERUM	RF
THYROID PEROXIDASE Abs	437	SST	SERUM	RF
THYROID PNL, SST	39051	SST	SST	
THYROID STIMULATING Ig	207L	SST	SERUM	RF
THYROTROPIN RECEPTOR Ab	4345L	RED SST	SERUM	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
THYROXINE BINDING GLOBULIN	4119L	RED SST	SERUM	RF
TIBC	157L	SST	SERUM	RF
TIBC CALCULATED	156	SST	SERUM	RF
TIBC CALCULATED	TIBC	SST	SERUM	RF
TIME COLLECTED	121L	SC	SEMEN	RM
TIME RECEIVED	122L	SC	SEMEN	RM
TNF-ALPHA	4520L	RED SST	SERUM	RF
TOBRAMYCIN, PEAK	027	RED	SERUM	RF
TOBRAMYCIN, PEAK	027L	RED	SERUM	RF
TOBRAMYCIN, RANDOM	030	RED	SERUM	RF
TOBRAMYCIN, RANDOM	030L	RED	SERUM	RF
TOBRAMYCIN, TROUGH	029	RED	SERUM	RF
TOBRAMYCIN, TROUGH	031L	RED	SERUM	RF
TOPIRAMATE, SERUM	1035L	RED	SERUM	RF
TORCH PROFILE, IgG	1009	RED,SST	SERUM	RF
TOTAL	TOTA	RED	SERUM	RF
TOTAL (IMI + DES)	2004L	RED,LAV,GREEN	SERUM/PLASMA	RF
TOTAL BILIRUBIN	127J	SST	SERUM	RF
TOTAL PA+NAPA	597L	RED	SERUM	RF
TOTAL PROTEIN	157	SST	SERUM	RF
TOTAL PROTEIN	157H	SST	SERUM	RF
TOTAL PROTEIN	6009B	SST	SERUM	RF
TOTAL PROTEIN, 24HR URINE	617	UC,24H U	URINE	RF
TOTAL PROTEIN, SERUM	157J	SST	SERUM	RF
TOTAL PROTEIN, URINE	616	SC	URINE	RF
TOTAL PROTEIN, URINE 24 HOUR	194	UC,24H U	URINE	RF
TOTAL VOLUME, URINE	512	UC,24H U	URINE	RF
TOXGRAN	TOXG	LAV	WHOLE BLOOD	RF
TOXOPLASMA Ab, IgG	851	SST	SERUM	RF
TOXOPLASMA, IgM	853L	SST	SERUM	RF
TRAMADOL (ULTRAM)	865M	SC	URINE	RF
TRANSFERRIN	112	SST	SERUM	RF
TRANSFORMATIONS PANEL	91401	SC	URINE	RF
TRANLYCYPROMIDE (PARNATE)	897M	RED GREEN	SERUM PLASMA	RF
TRAZODONE (DESYREL)	866M	SC	URINE	RF
TREASURE COAST DOT DRUG SCREEN	2040	SC	URINE	RF
TREATMENT UDS(13)	70770	SC	URINE	
TRICHO	U24	SC	URINE	RM
TRICHOMONAS CULTURE	118L	FLUID	SWAB	RM
TRICHOMONAS CULTURE	118Q	"INPOUCH"	FLUID	RF
TRICHOMONAS SCREEN	980	FLUID	SWAB	RM
TRICYCLIC ANTIDEPRESSANTS	411L	SC	URINE	RF
TRICYCLIC CONF (GC/MS)	571QU	SC	URINE	RF
TRICYCLICS	571	SC	URINE	RF
TRICYCLICS	571Q	SC	URINE	RF
TRIGLYCERIDE	158H	SST	SERUM	RF
TRIGLYCERIDES	158	SST	SERUM	RF
TRIGLYCERIDES	146L	SST	SERUM	RF
TRIGLYCERIDES	4041L	SST	SERUM	RF
TROPONIN I	381	SST	SERUM	RF
TROPONIN I	381B	SST	SERUM	RF
TROPONIN I	381L	SST	SERUM	RF
TROVATO CTR MALE PNL-3SST,LAV,R	91291	SST,L,R	BLD, PLASMA	
TROVATO THYROID PNL-2SST	91290	2SST	SERUM,	
TRTMNT FEMALE ADM PNL-2SST,LAV	70772	2SST,LAV	SERUM, BLD	
TRTMNT MALE ADM PNL-2SST,LAV	70771	2SST,LAV	SERUM, BLD	
TRYPSIN	702L	SST	SERUM	FZ
TRYPTASE	2021L	RED SST	SERUM	RF
TSH	113	SST	SERUM	RF
TSH	4204L	SST	SERUM	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
t-TRANSGLUTAMINASE IgA	643L	SST	SERUM	RF
t-TRANSGLUTAMINASE IgG	677L	SST	SERUM	RF
t-TRANSGLUTAMINASE, IgA	676L	SST	SERUM	RF
uE3 Mom	461L	SST	SERUM	RF
uE3 VALUE	459L	SST	SERUM	RF
UIBC	156L	SST	SERUM	RF
UNSTIMULATED CD4	1096L	GREEN	WHOLE BLOOD	RT
UR-AC	U26	SC	URINE	RF
UREA NITROGEN, 24HR URINE	926L	24 HR U.C.	24 HR URINE	RF
UREA NITROGEN, URINE	925L	SC	URINE	RF
UREAPLASMA/MYCOPLASMA CULTURE	4347L	VTM	FLUID	RF
URIC ACID	159	SST	SERUM	RF
URIC ACID	159H	RED SST	SERUM	RF
URIC ACID, 24HR URINE	928L	24 HR UC	24 URINE	RF
URIC ACID, URINE	927L	SC	URINE	RF
URINALYSIS	040	SC	URINE	RF
URINALYSIS	040H	SC	URINE	RF
URINALYSIS	040MM	SC	URINE	RF
URINALYSIS W/RFX URINE CULTURE	043	SC	URINE	RF
URINALYSIS WITH MICROSCOPIC	056	SC	URINE	RF
URINE ADMIT	2012	SC	URINE	RF
URINE ADMIT PANEL	1061	SC	URINE	RF
URINE CULTURE	090	SC	URINE	RF
URINE CYTOLOGY	473L	SC	URINE	RF
URINE CYTOLOGY	473Q	SC	URINE	RF
URINE GLUCOSE 1 HR	UGL1H	SC	URINE	RF
URINE GLUCOSE 1/2 HR	UGL15	SC	URINE	RF
URINE GLUCOSE 2 HR	UGL2H	SC	URINE	RF
URINE GLUCOSE 3 HR	UGL3H	SC	URINE	RF
URINE GLUCOSE 4 HR	UGL4H	SC	URINE	RF
URINE GLUCOSE 5 HR	UGL5H	SC	URINE	RF
URINE GLUCOSE 6 HR	UGL6H	SC	URINE	RF
URINE PROTEIN ELECTROPHORESIS, URINE	112Q	SC RANDOM	URINE	RF
URINE SATURATION	4659L	24HR UC	URINE	RF
UROBILINOGEN	U10	SC	URINE	RF
UROVYSION FISH, BLADDER CANCER	4326L	TCC MONITOR KIT	URINE	RF
VACUOLE	VACU	LAV	WHOLE BLOOD	RF
VAGINAL CULTURE	086	CULTURETTE,SWAB	SWAB	RM
VAGINITIS/VAGINOSIS DNA PROBE	4295L	SWAB VP III	SWAB VP III	RT RF
VALIDITY TESTING (TEST)	99998			
VALIDITY TESTING PANEL	539	UC	URINE	RF
VALIDITY TESTING PANEL (TEST)	9998			
VALPROIC ACID	079	RED	SERUM	RF
VALPROIC ACID	079MM	RED	SERUM	RF
VANCOMYCIN PEAK	038A	RED	SERUM	RF
VANCOMYCIN RANDOM	038C	RED	SERUM	RF
VANCOMYCIN TROUGH	038	RED	SERUM	RF
VANCOMYCIN, PEAK	038AMM	RED	SERUM	RF
VANCOMYCIN, TROUGH	038MM	RED	SERUM	RF
VAP CHOLESTEROL PANEL	4086L	SST	SERUM	RF
VARICELLA ZOSTER Ab, IgM	592L	SST	SERUM	RF
VARICELLA ZOSTER IgG	402	SST	SERUM	RF
VENIPUNCTURE	001			
VENIPUNCTURE	1			
VENIPUNCTURE/STAT	002			
VENLAFAXINE (EFFEXOR)	923M	RED GR	SERUM PLASMA	RF
VGCC ANTIBODY	4212L	RED SST TO TT	SERUM	RF
VIP, PLASMA	5236L	LAV	PLASMA	FZ
VIRAL CULT, RAPID RFX H1N1	4586L	SWAB VTM	SWAB	RF
VIRAL CULTURE	5043Q	SC	FLUID TISSUE	RT

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
VIRAL CULTURE, GENERAL	5043L	SWAB VTM	SWAB VTM	RF
VIRAL CULTURE, RAPID INFLUENZA	159L	VTM	FLUID	RM
VIRAL CULTURE, RAPID, LESION	4597L	VTM	FLUID SWAB	RF
VISCOSITY, SERUM	4533L	RED SST	SERUM	RF
VITAMIN A, SERUM	385L	RED SST FOIL	SERUM	RF PROTECT FROM LIGHT
VITAMIN B1	378L	LAV	PLASMA	FR
VITAMIN B12 UNSAT BIND CAP	4709L	RED SST LAV	SERUM PLASMA	RF
VITAMIN B2	4686L	LAV FZ, FOIL	WHOLE BLOOD	FZ
VITAMIN B6	290L	LAV	PLASMA	FZ
VITAMIN C	4561L	GR RED SST FOIL	PLASMA/SERUM	FZ FOIL
VITAMIN D 1,25 DIHYDROXY	556Q	RED SST TO TT	SERUM	RF
VITAMIN D, 1,25 DIHYDROXY	556L	SST	SERUM	RF
VITAMIN D, 25-HYDROXY	460	RED SST	SERUM	RF
VITAMIN E, SERUM	4302L	RED SST	SERUM	RF(PROTECT FROM LIGHT)
VITAMIN K1, SERUM	868L	SST	SERUM	FZ
VLDL (calc)	4043L	SST	SERUM	RF
VLDL (CALCULATED)	345	SST	SERUM	RF
VLDL calc	4077L	SST	SERUM	RF
VMA, 24 HR URINE	579L	24 HOUR URINE	URINE	RM
VMA, URINE RANDOM	578L	SC	URINE	RM
VOLUME	125L	SC	SEMEN	RM
vWF ACTIVITY	191L	BLUE	PLASMAN	FZ
vWF ANTIGEN	190L	BLUE	PLASMA	FZ
vWF PROFILE (VIII)	188L	BLUE	PLASMA	FZ
W001 RAGWEED, SHORT / COMMON	2075L	SST	SERUM	RF
W005 WORMWEED	3021L	SST	SERUM	RF
W009 PLANTAIN, ENGLISH	2076L	SST	SERUM	RF
W010 LAMB'S QUARTER	2077L	SST	SERUM	RF
W014 PIGWEED, ROUGH	3022L	SST	SERUM	RF
W015 LENSACLE	4682L	SST	SERUM	RF
W018 SHEEP SORREL (DOCK)	2078L	SST	SERUM	RF
W020 NETTLE	3023L	SST	SERUM	RF
W023 DOCKWEED, YELLOW	3024L	SST	SERUM	RF
W028 ROSE	4336L	RED SST	SERUM	RF
W046 FENNEL, DOG	2079L	SST	SERUM	RF
W067 BACCHARIS (GROUNDSEL)	2080L	SST	SERUM	RF
W069 IODINE BUSH	4213L	SST	SERUM	RF
W075 WING SCALE	3025L	SST	SERUM	RF
WALKER SMAC PANEL	3004	2SST	SERUM	RM
WARFARIN, SERUM	4590L	RED LAV GR	SERUM PLASMA	RF
WATERSHED ADM FEMALE (3SST LV)	71242			
WATERSHED ADM MALE(3SST LAV)	71241			
WATERSHED UDS PANEL (UA)	71240	SC	URINE	RF
WBC	8000B	LAV	WHOLE BLOOD	RF
WBC	WBC	LAV	WHOLE BLOOD	RF
WBC	WBCJ	LAV	WB	RF
WBC'S	U13	SC	URINE	RM
WEIGHT	452L	SST	SERUM	RF
WEISS PANEL II	1033	SST,L	SERUM/PLASMA	RM
WEISSMAN PANEL 1	2019	SST,L	SERUM/PLASMA	RM
WEISSMAN PANEL 2	3000	SST,L	SERUM/PLASMA	RM
WEISSMAN RHEV.PANEL	3003	2SST	SERUM	RM
WEST NILE VIRUS ANTIBODY, SERUM	983L	RED SST	SERUM	RF
WEST NILE VIRUS, IgG	981L	RED SST	SERUM	RF
WEST NILE VIRUS, IgM	982L	RED SST	SERUM	RF
WESTERN BLOT	5223L	SST	SERUM	RF
WET PREP	5044L	INPOUCH TV	SWAB FLUID	RT
WET PREP	HC3	INPOUCH TV	SWAB FLUID	RT
WOUND CULTURE	612	SWAB	FLUID	RM
XYLOSE TOLERANCE	4350L	GRAY	PLASMA	RF

DESCRIPTION	CODE	CONTAINER	SPECIMEN	TEMPERATURE
YEAST	U20	SC	FLUID	RM
ZIELINSKI PANEL	2075	SST,R,2L	SERUM/PLASMA	RM
ZINC PROTOPORPHYRIN	151L	LAV/RB	WHOLE BLOOD	RM
ZINC, BLOOD	386L	RED OR RB	SERUM/PLASMA	RT
ZINC, URINE	4528L	UC OR 24HRUC	URINE	RF
ZOLPIDEM (AMBIEN)	941M	RED GREEN	SERUM PLASMA	RF
ZONISAMIDE (ZONEGRAN)	291L	RED OR LAB	SERUM/PLASMA	RF